Keralloy® KB

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Keralloy® KB

Instruction for use

Area of application: Non precious casting alloy based on cobalt for metal-ceramic accord. DIN EN ISO 22674, type 5.

The product is intended for dental laboratories for the creation of products that must be rigid and resistant, single crowns, bridges and superstructures on implants. The alloy Keralloy® KB is distinguished by an outstanding fluidity, which grants to fill even the thinnest details of the frame, down to three tenths of millimeter. Their molecular structure allows you to obtain smooth, compact surfaces with little oxide formation even after repeated oven treatment and their thermal expansion coefficient is ideal for every ceramic of the latest generation. The alloy Keralloy® KB is produced in conformity with standards ISO 9693:2020 and ISO 22674:2016, it is highly corrosion resistant, it does not contain Nickel and it is completely free of toxic elements like beryllium, indium and gallium. Keralloy® KB is a product intended to be used by qualified technician.

Intended purpose: Fabrication of metal base structures for metal-ceramic crowns and FPDs.

Intended patient population: This alloy is intended to be used for oral rehabilitation in edentulous or partially edentulous patients who suffer from loss of hard dental tissue. There are no limitations related to age, sex, clinical condition or other conditions, other than those possibly determined case by case by the professional user.

Intended users: The product is intended for professionals of the dental sector, in particular: prosthodontic laboratories, dental technicians, dentists. These professional figures shall have a high degree of specialisation and are perfectly aware of the destination of the product and the correct method of use for the realization of a final product and elaborate different methods to reach the best shape of the final prosthesis. Contact and/or manipulation by other people is not envisaged.

Instruction for use

Modeling: Modeling with a minimum thickness of 0.3mm on simple structures and 0.5mm in complex cases or with patients with bruxism. Avoid elbows.

Casting Sprues: Direct pinning: Pin: Ø 3mm, Nourice: Ø >= 6mm

Indirect pinning with the stabilizer bar: Pin: Ø 3mm, Bar >= 5mm, Casting sprue >= 6mm

Casting: Cast Keralloy® KB alloys in pre-heated and cleaned ceramic crucibles used for this alloy only. Each alloy should have its crucible; do not overheat the alloy and do not use the flux. Reach a temperature of 850-950°C; maintenance time in pre-heated oven depends on investment characteristics and on the dimension of the cylinder. Induction casting: start the casting when the cylinders will melt together and just before the superficial crust opens. Open flame melting: place ingots in the pre-heated ceramic crucible and heat them evenly with circular movements. When ingots have melted, start centrifugal unit. Use multi-flame welding torches only. Do not use any flux. Indicative values for flame regulation: Acetylene 0,4 bar / Oxygen 2 bar; Propane 0,2 bar / Oxygen 2 bar.

Manufacturing: Let the cylinder cool down at room temperature, remove the investment and sandblast with aluminium oxide 110 to 150 µm at a pressure of 3-4 bar. It is recommended an homogenization and degassing treatment of the casted piece at 1000°C for 10 minutes in atmosphere. Continue manufacturing process using tungsten carbide burs. The worked framework has to be sandblasted with disposable equipment in aluminium oxide of 110 to 150 µm at a max. pressure of 3-4 bar and then steam clean. After cleaning, the framework should not be touched anymore and should be held with clamp forceps only.

Oxidation: Fire for 5 min. under vacuum at 950-980°C. After firing, the oxide layer has to be carefully sandblasted using disposable oxide-aluminium equipment 110 to 150 µm at a pressure of 2.5-4 bar. Steam and check that the surface has an homogeneous grey surface.

Re-use: The device can be reused until the available material is exhausted.

Re-use of sprues: The best results are obtained with pure Keralloy® KB. It is nonetheless possible to re-use sprues and cones only once, given that new metal coming from the same lot is added in equal quantity.

Ceramization: It's possible to use Bonding. Apply the ceramic, following its manufacturer's instructions. Slow cooling is recommended.

Soldering: We suggest using SILADENT Cr-Co solder to weld.

Waste disposal: The processing scrap must be disposed of as special waste in accordance with the EC directives 2008/98/CEE on waste, and 94/62/CEE on packaging and packaging waste and in compliance with national legislation in force on the subject.

Safety instructions

- Cobalt-Chrome-based alloys can seldom cause dermatitis on sensitive subjects. A Patch-Test is thus advisable.
- Metal dusts and smoke are dangerous for health. Use exhaust fans while casting and polishing.
- The alloy contains Cobalt, element classified as CMR (CAS Number: 7440-48-4).
- Before prosthesis application verify if other metal elements are in patient's oral cavity. Coexistence of different metals can cause a "pile" effect.
- We identify every batch with a number. We recommend to write it down in patient's file to allow its complete traceability.
- This product does not need any special preservation precautions. MESA keeps batch traceability for 15 years.
- The methods of preserving the products made with Keralloy® KB are under the responsibility of professional operators (dental technicians). This product must be stored, transported and used at a temperature between -30°C and +50°C to ensure the integrity of the packaging and the instructions contained therein.
- MESA recommends that the patient be made aware of the possibility for dental alloys to affect MRI results.
- The product is sold non-sterile.

Users shall report any major accident occurred in relation to the use of this medical device to the manufacturer and the competent authority of the Member State in which the user and/or patient is established.

Composition: (in % by mass)		Technical data: (guidelines)	
Co	64,0	Proof stress 0,2% (MPa)	570
Cr	21,0	Elongation at rupture (%)	10
Mo	6,0	Modulus of elasticity (GPa)	194
W	6,0	Coefficient of expansion	
		25-500 °C	14,1x10 ⁻⁶ K ⁻¹
		25-600 °C	14,6x10 ⁻⁶ K ⁻¹
Other constituents: Si, Fe, Mn		Tensile strength (MPa)	734
		Vickers hardness VH 10	286
		Density (g/cm ³)	8,8
Tolerances in the composition of the alloy (% by mass) are within the range of current DIN standards.		Melting range (°C)	1.309-1.417
		Casting temperature (°C)	1.460
		Max. oxide firing temperature (°C)	935

Symbol	Description	Symbol	Description
	Batch number		Date of manufacture
	Consult instructions for use		
	Manufacturer		CE Marking
	Weight estimate		Non-sterile
	Medical device		Catalogue number
	Temperature limitation		Contains hazardous substances
	Product UDI Code		

Mesa Italia S.r.L. via dell'Artigianato, 35/37, 25039 Travagliato (BS)
ITALY Tel. +39 0306863251
info@mesaitalia.it www.mesaitalia.it

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SILADENT Dr. Böhme & Schöps GmbH