1. Intended use / Indication

SprintRay EU High Impact Denture Base is a light-curable, polymerizable resin to be used in conjunction with extraoral light-curing equipment.

SprintRay EU High Impact Denture Base is indicated for the fabrication and repair of removable full and partial dentures and base plates by additive manufacturing.

2. Contraindication

SprintRay EU High Impact Denture Base restorations are contraindicated ...

1. ... if a patient is known to be allergic to any of the ingredients.
2. ... direct intraoral contact of uncured or partially uncured material.
3. ... for every application that is not part of the indication (see above).

3. Patient Target Group

Persons being treated in the context of a dental procedure.

4. Intended Users

Dentist, dental technician

5. Requirements

Software:

- RayWare (CAD)

Hardware:

- Pro Desktop 3D (3D Printing)

- Pro Wash / Dry (Cleaning)

- ProCure (Post Curing)

SprintRay, Inc., 3577 North Figueroa St., Los Angeles, CA 90065

6. Material

SprintRay EU High Impact Denture Base consists of functional (meth)acrylic resins, initiators, dyes and stabilisers.

7. Geometric presetting

Lingual Ridge ≥ 2.5 mm (Maxilla / Mandibula)

Palatal/Lingual ≥ 3.0 mm (Maxilla) / ≥ 2.5 (Mandibula)

Facial/Buccal ≥ 2.5 mm (Maxilla / Mandibula)

8. Material parameters

Radiation penetration depth controllable by exposure time

≥ 50 µm (4 mils)

≥ 100 µm (5 mils)

9. Manufacturing process (fig. 1-10)

1. Prepare data (RayWare).
2. Choose process parameters (build style, etc.).
3. Transfer prepared data to 3D printing device.
4. Prepare 3D printing – shake bottle.
5. Fill resin tank of 3D printing device.
6. Build the parts.
7. Clean parts (ProWash / Dry): 5 min. / 5 min.
8. Dry parts (ProWash / Dry): 3 min.
9. Post curing (ProCure): 5 min. / 60 °C.
10. Finish parts.

10. Bonding of denture teeth


If using printed denture teeth – see additional IFU **SprintRay EU Temporary Crown & Teeth** for details. In order to realize the optimal bonding of artificial teeth to the dental object, it can be necessary to roughen the surface of the denture and to use an appropriate adhesive. The suitability of a bonding agent or adhesive must be approved prior use. Recommended bonding agent is VITA VIONIC® BOND. The commonly used techniques and dental work steps must be carefully considered when manufacturing the dentures. Please follow the instructions and indications of the corresponding suppliers.

11. Notice

Please follow the instructions provided by the software manufacturer in terms of parameter settings and design recommendations. Please follow the instructions provided by the hardware manufacturer in terms of parameter settings/printing & post curing recommendations. To avoid detrimental effects on material quality do not expose the liquid material to irradiation under any circumstances. Deviations from the described manufacturing processes or storing conditions may lead to different mechanical and optical properties of the material. Ensure personal protective gear during processing. According to the EU Medical Devices Regulation, users/patients are obliged to report serious events with a medical device to the manufacturer and to the competent authority of the country in which they occurred. **Caution:** The lot number and the best before date are indicated on each packaging. In case of claims please always indicate the lot number of the product. Do not use the product after expiry of the best before date. Avoid breathing dust/fume/gas/mist/vapours/spray. Wash hands thoroughly after handling. Use only outdoors or in a well-ventilated area. Contaminated work clothing should not be allowed out of the workplace. Wear protective gloves/protective clothing/eye protection/face protection. IF ON SKIN: Wash with plenty of soap and water. IF INHALED: Remove person to fresh air and keep comfortable for breathing. IF IN EYES: Rinse cautiously with water for several minutes. Remove contact lenses, if present and easy to do. Continue rinsing. Call a POISON CENTER or doctor/physician if you feel unwell. Specific treatment (see information on this label). If skin irritation or rash occurs: Get medical advice/attention. If eye irritation persists: Get medical advice/attention. Take off contaminated clothing and wash before reuse. Store in a well-ventilated place. Keep container tightly closed. Store locked up. Dispose of contents/container according to official regulations.

12. Hazard statements

Causes skin irritation. May cause an allergic skin reaction. Causes serious eye irritation. May cause respiratory irritation.

Physical properties*/
Physikalische Eigenschaften*/
Propriétés physiques*/
Propiedades físicas*/
Proprietà fisiche*:

SprintRay EU High Impact Denture Base

■ Shore hardness D/
Shore-Härte D/
Dureté Shore D/
Durezza Shore D/
Durezza Shore D
(ISO 48-4):
> 80**

■ Flexural strength/
Biegefestigkeit/
Résistance à la flexion/
Fuerza flexible/
Resistenza alla flessione
MPa (ISO 20795-1):
≥ 65

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SprintRay EU High Impact Denture Base

1 kg:
λ ≤ 405 nm

orange-pink/
orange-pink/
orange-rose/
rosa-naranja/
rosa-arancione
REF: SRE-1001311

light-pink/
hell-pink/
rose-clair/
rosa-claro/
rosa-chiaro
REF: SRE-1001312

deep-pink/
dunkel-pink/
rose-profond/
rosa-profundo/
rosa-scuuro
REF: SRE-1001313

light-pink B/
hell-pink B/
rose-clair B/
rosa-claro B/
rosa-chiaro B
REF: SRE-1001314

natural-pink/
natur-pink/
rose-naturel/
rosa-natural/
rosa-naturale
REF: SRE-1001315

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5. Anforderungen
Software:
- RayWare (CAD)

Hardware:
- Pro Desktop 3D (3D Printing)
- Pro Wash / Dry (Cleaning)
- ProCure (Post Curing)

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