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Stand der Information: 05/2017

Druckdatum: 15.05.2017

CoNiCr Special solder Instructions for use

CE 0124

Area of application: Non precious cobalt based hard solder alloy acc. EN ISO 9333.

Product description: SILADENT CoNiCr special solder is indicated to weld cobalt based alloys.

REF 102877 CoNiCr Special solder, 1 mm (10 g)

This special solder contains 21 % nickel and should not be processed by individuals with known nickel sensitivity (allergy). Metal dust contains Ni and must not be inhaled.

Working instruction: Objects which have to be soldered must be clean and free of all oxides. Then cover the area to be soldered with flux. Mix the flux with water to a pasty consistence. Finally heat up the area slowly with a torch and deposit the solder evenly within the soldering gap.

Composition (in % by mass): Co: 52 Ni: 21 Cr: 20

Other constituents: Fe-Si-C-Mn

Melting range (°C): Solidus-/liquidus temperature: 992 – 1.185 °C

Flow temperature: 1.240 °C

Safety hint: Metal dust is harmful to your health. When deflasking and blasting use a suction extraction system and breathing mask!

Contraindications, side effects and interactions with other dental alloys: In cases of hypersensitivity (allergy) to the constituents of the alloy. Hypersensitivity reactions (allergies) or electrochemically-induced local dysaesthesia have been reported in individual cases. A Patch-Test is thus advisable. Galvanic effects can occur under proximal or antagonistic contact with dentures of different alloys. The alloy can affect MRI results. Each supply is identified by a lot-number. Take note of this number on the patient's file in order to trace down the product.

Storage conditions: No special storage conditions are required.

Further questions: Contact our technical services (Tel.: +49(0)53 21–37 79 25/26) or our sales representatives.

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