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**Area of application:** Non precious casting alloy based on cobalt for metal-ceramic accord. DIN EN ISO 22674, type 4.

Composition: (in % by mass)		Technical data: (guidelines)	
Co	63,4	Proof stress 0,2% (MPa)	577
Cr	28,85	Elongation at rupture (%)	16
Mo	6,1	Modulus of elasticity (GPa)	210
<b>Other constituents:</b> Si, Mn		Coefficient of expansion 25-500 °C	14,7x10 <sup>-6</sup> K <sup>-1</sup>
Tolerances in the composition of the alloy (% by mass) are within the range of current DIN standards.		Tensile strength (MPa)	830
		Vickers hardness VH 10	310
		Density (g/cm <sup>3</sup> )	8,4
		Melting range (°C)	1.370-1.430
		Casting temperature (°C)	1.500
		Max. oxide firing temperature (°C)	1.050

### Investing and casting:

Keralloy® FG can be used with phosphate-bonded crown and bridge investments, such as Premium, Presto Vest II and TeleVest. It is preheated to 850 °C and heat soaked for at least 30 minutes before casting. Follow the recommendations in the SILADENT technique instructions when attaching sprues. Use a separate ceramic crucible for casting Keralloy® FG. High frequency casting centrifuge: The alloy is cast when all the ingots have melted together and a final shadow runs over the molten metal just before the oxide layer disintegrates. Please note that the melting point may vary depending on the type of casting machine used. During the flame casting wait after the ingots collapse until the solidified melt point move spherical in the crucible, then cast immediately.

### Prepare and polishing:

After sandblasting, prepare the metal surface with corundum, diamond or carbide milling tools. For polishing of Keralloy® FG don't use fine corns sized rubbers (wheels or roller). Afterwards polish with goats hair brushes the surfaces, chewing surface should be polished with metal brushes under use of a suitable polishing paste. You will reach a gloss finish with linen or cotton buff under usage of high gloss polishing paste or similar polishing materials.

### Firing porcelain:

Sandblast the metal surface with 250 µm aluminium oxide. It's recommended the pickle the metal with ethyl acetate for approx. 3 minutes in a ultrasonic and afterwards steam it. It is important that the complete frame must be covered thru the ethyl acetate. After that start with the oxide firing of the porcelain. All oxide firings (except the opaque firing) must be handled with a compensation cooling. Clean the frames after every firing under rising water (brush it off).

### Soldering:

Clean the metal surface without sandblasting. Heat up the soldered object to the temperature of the solder (e.G. SILADENT CoCr solder). If the necessary temperature is reached, the solder should flow into the soldered joint.

### Safety hint:

Metal dust is harmful to your health. When preparing and blasting use a suction extraction system and breathing mask!

### Contraindications, side effects and interactions with other dental alloys:

In cases of hypersensitivity (allergy) to the constituents of the alloy, discontinue its use. In individual cases, hypersensitivity reactions (allergies) and electrochemically induced local dysaesthesia have been reported, such as changes in taste and irritation of the oral mucosa. Galvanic effects can occur under proximal or antagonistic contact with dentures of different alloys.

### Storage conditions:

No special storage conditions are required.

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