# Keralloy® KB

## CE 0124

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## CE 0124

### Instructions for use

**Area of application:** Non precious casting alloy based on cobalt for metal-ceramic accord. DIN EN ISO 22674, type 5.

| Composition: (in % by mass)   |      | Technical data: (guidelines)       |                                       |
|---|------|------------------------------------|---------------------------------------|
| Co  | 64,0 | Proof stress 0,2% (MPa)            | 570                                   |
| Cr  | 21,0 | Elongation at rupture (%)          | 10                                    |
| Mo  | 6,0  | Modulus of elasticity (GPa)        | 194                                   |
| W   | 6,0  | Coefficient of expansion           |                                       |
|   |      | 25-500 °C                          | 14,1x10 <sup>-6</sup> K <sup>-1</sup> |
|   |      | 25-600 °C                          | 14,6x10 <sup>-6</sup> K <sup>-1</sup> |
| <b>Other constituents:</b> Si, Fe, Mn   |      | Tensile strength (MPa)             | 734                                   |
|   |      | Vickers hardness VH 10             | 286                                   |
|   |      | Density (g/cm <sup>3</sup> )       | 8,8                                   |
| Tolerances in the composition of the alloy (% by mass) are within the range of current DIN standards. |      | Melting range (°C)                 | 1.309-1.417                           |
|   |      | Casting temperature (°C)           | 1.460                                 |
|   |      | Max. oxide firing temperature (°C) | 935                                   |

#### Product description:

Keralloy® KB is a non-precious, cobalt based bonding alloy. Cause of the technical properties very thin-walled constructions are possible. The structure of the alloy generates less oxidation also after repeated casting and offers a gut processing and polishing. Its thermal expansion coefficient is ideal for every ceramic of the latest generation. Keralloy® KB is highly corrosion resistant and does not contain any beryllium, indium or gallium ( $\leq 0,1\%$ ).

#### Investing and casting:

Keralloy® KB can be used with phosphate-bonded crown and bridge investments, such as Premium, Presto Vest II and TeleVest. It is preheated to 850 °C and heat soaked for at least 30 minutes before casting. Follow the recommendations in the SILADENT technique instructions when attaching sprues. Use a separate unglazed and preheated ceramic crucible for casting Keralloy® KB to prevent it being contaminated by other alloys. Clean the crucible after each cast. With the high-frequency melting the alloy is cast when all the cylinders have melted together and a final shadow runs over the molten metal just before the oxide layer disintegrates. Please note that the melting point may vary depending on the type of casting machine used and when melting with vacuum, for example, less oxide forms and the oxide layer in the vacuum disintegrates more quickly. In general cast with incorporating flux. Considerable experience is required when melting high-grade alloys and it is most important that heating is set exactly according to the manufacturer's instructions. Sample values for flame regulation: Acetylene 0,4 bar / Oxygen 2 bar; Propane 0,2 bar / Oxygen 2 bar; Leading pressure methane / Oxygen 2 bar. Inaccurate heat setting impairs the alloy and can have a detrimental effect on the porcelain bond. Place ingots in the pre-heated ceramic crucible and heat them evenly with circular movements. When ingots have melted, start centrifugal unit. Use multiflame welding torches only. Do not use any flux. Casting Keralloy® KB repeatedly is not recommended, as this could alter both its composition and properties. It is nonetheless possible to re-use sprues and cones only once, given that new metal coming from the same lot is added in equal quantity.

#### Manufacturing and polishing:

Cool down the muffle at room temperature, remove the investment and sandblast with aluminum oxide 100 µm to 250 µm at a pressure of 3-4 bar. Continue manufacturing process using fraises for hard metal.

#### Firing porcelain:

1. The worked framework has to be sandblasted with disposable equipment in aluminum oxide of 110 to 150 µm at a max. pressure of 3-4 bar and then steamclean. After cleaning, the framework should not be touched anymore by hand.
2. Fire for 5 min. under vacuum at 950 - 980 °C (10 °C more than back-ground material firing temperature).
3. After firing, the oxide layer has to be carefully sandblasted using disposable oxide-aluminum equipment 110 to 150 µm at a pressure of 2,5 - 3 bar, then steam.
4. If the framework has a homogeneous grey surface, then you can apply and fire ceramic, following its manufacturer's instructions.
5. It is possible to use all ceramics normally available on the market for CoCr alloys. Apply manufacturing instructions of the ceramic producer (e.g. VITA VM13). Do not use any bond.

#### Soldering and welding:

Presolder using SILADENT Cr-Co solder and the appropriate high-fusing flux or a high-fusing solder for precious bonding alloys. Laser weld using a SILADENT Cr-Co laser welding rod (REF 102806).

#### Safety hints:

Metal dust is harmful to your health. When preparing and blasting use a suction extraction system and breathing mask!

#### Contraindications, side effects and interactions with other dental alloys:

In cases of hypersensitivity (allergy) to the constituents of the alloy, discontinue its use. A Patch-Test is thus advisable. In individual cases, hypersensitivity reactions (allergies) and electrochemically induced local dysesthesia have been reported, such as changes in taste and irritation of the oral mucosa. Galvanic effects can occur under proximal or antagonistic contact with dentures of different alloys. It is recommended that the patient be made aware of the possibility for dental alloys to affect MRI results. Each supply is identified by a lot-number. Take note of this number on the patient's file in order to trace down the product.

#### Storage conditions:

No special storage conditions are required.

**Further questions:** Contact SILADENT technical services  
(Tel.: +49 (0) 53 21 - 37 79 25 / 26)  
or our sales representatives.

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