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Modiral® S Instructions for use

Area of application: Non precious alloy based on cobalt for removable dentures acc. DIN EN ISO 22674, type 5.

Composition: (in % by mass)		Technical data: (guidelines)	
Co	62,5	Density ρ (g/cm ³)	8,2
Cr	30,0	Vickers hardness HV 10	375
Mo	5,1	Melting range $T_S - T_L$ (°C)	1.260 - 1.320
Si	1,0	Casting temperature T_{cast} (°C)	1.440
Mn	1,0	0,2-% Yield strength $R_{p,0.2}$ (MPa)	620
Other constituents: (< 1%) C, N, Nb, Fe		Modulus of elasticity (GPa)	220
Tolerances in the composition of the alloy (% by mass) are within the range of current DIN standards.		Tensile elongation at break A_5 (%)	4,5

Product description:

Modiral® S is an extra hard partial denture alloy based on cobalt and finds application in frame works. Modiral® S is free from nickel, cadmium, beryllium and lead and fulfils the standards of EN ISO 22674 type 5 for appliances in which parts require the combination of high stiffness and strength, e.g. thin removable partial dentures, parts with thin cross-sections and clasps.

Waxing-up:

Lead round wax wires with Ø 3,5-4 mm in direction of flow, avoid sharp deflections.

Investing and casting:

Suitable investment materials are phosphate-bounded investments for partial denture alloys. Preheat the investment to 950-1000 °C. Depending on the model and the casting machine and when total plates are modelled the final temperature is 1050 °C. Hold final temperature for 45-60 minutes depending on the size of the investment ring and degree of filling. Follow the manufacturer's instructions for use of the casting machine. Always use an individual ceramic crucible for Modiral® S to prevent contamination with other alloys. Clean crucible after each use to avoid residues of slag. Do not overheat the alloy. Start casting as soon as the ingots have collapsed giving a uniform melt. For melting by flame heat the ingots and give a rotary motion by use of the flame. Start casting as soon as the bath begins to vibrate. Allow the cylinder to cool down slowly to the ambient temperature and deflask without hitting the cone.

Treatment of the Cast Object:

After deflasking:

1. Sand blast the surface by use of a pencil-blaster with aluminium oxide 100 µm or 250 µm.
2. Separation of sprues and finishing of the object.
3. Electrolytic polishing with commercial electrolytes in dental polishing units. Cover clamps and fitting parts with covering varnish while polishing.

Finishing

After finishing and fitting smooth the frame with a rubber polisher.

Soldering and welding:

Soldering can be carried out with commercially available solders und high temperature flux. The width of the solder gap should be 0,05-0,2 mm. For welding with laser use suitable commercially available metal welding wires.

Safety note:

Metal dusts are harmful to health. Use a dust extractor. Consider allergic hypersensitivities to contents of the alloy. In case of suspected incompatibility with individual elements of this alloy, this should not be used.

Other notes:

All serious incidents related to the device shall be reported to the manufacturer and to the competent authority of the member state in which the user and/or the patient is established.

Warranty:

These application recommendations are based on own experiments and experiences and can therefore only be regarded as guidelines. The dentist or dental technician is responsible for the correct processing of this alloy.

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Manufacturer: German Special Alloys GmbH | Carl-Friedrich-Benz-Str. 1b | 47877 Willich | Germany

German Special Alloys packaging carries the mark **CE 0044**

Symbol	Description	Symbol	Description
	Batch number		Manufacturer
	Refer to instructions for use		Not for reuse

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Distribuidor: SILADENT Dr. Böhme & Schöps GmbH, Im Klei 26, DE 38644 Goslar
Fabricante: German Special Alloys GmbH | Carl-Friedrich-Benz-Str. 1b | 47877 Willich | Germany
Distributore: SILADENT Dr. Böhme & Schöps GmbH, Im Klei 26, DE 38644 Goslar
Fabbricante: German Special Alloys GmbH | Carl-Friedrich-Benz-Str. 1b | 47877 Willich | Germany