Presto Vest II

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DIN EN ISO 15912, Typ 1, Klasse 2

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DIN EN ISO 15912, type 1, class 2

Description:

Presto Vest II is a phosphate-bonded, special speed investment for crown and bridge work. A new form of grain distribution allows Presto Vest II to be placed in a preheated furnace at 850°C, just 15 minutes after investing. It is ready to cast after 45 - 60 minutes depending on the size of the ring.

Presto Vest II is graphite-free and suitable for casting all types of dental alloy.

Waxing up:

Optimum fit and surface quality are best achieved using patterns fabricated entirely from wax. Using resins often results in much poorer surface quality and cracks

in the mould with the speed cast technique. This applies particularly to waxing up removable sections of attachment and telescopic crown work.

Crown and bridge work:

Mix Presto Vest II according to the instructions and pour it up to the patterns without vibration. Then vibrate gently into the patterns using a small instrument, if necessary. Fill the ring completely without any further vibration. Adhere strictly to the time given in the instructions for the speed technique.

Use of pressure:

Investing under pressure is not necessary due to Premium's extremely homogenous grain distribution. If investing under pressure is preferred, the pressure should be released after 15 minutes so that it does not impede the setting expansion.

Preheating for speed casting:

Place the mould in a preheated furnace at 850°C, 15 minutes after mixing the investment. Ensure that there is a space of approx. 5.0 mm between the mould and the furnace floor by placing the mould on pieces of investment or stilts. If a different preheating temperature is required for the alloy being cast, the temperature can be adjusted after the first 15 minutes. The furnace fan and extractor should remain switched off during the first 15 minutes and the furnace door should not be opened because of the risk of instant combustion. The preheating time for size 3 rings is approx. 40 minutes and for size 6 rings approx. 60 minutes.

Technical Data:

Mixing liquid:	SILADENT type 100 expansion liquid
Mixing ratio Powder : liquid:	100 g : 24 ml
Manual mixing:	10 sec.
Mixing under vacuum:	60 sec.
Working time:	5 - 7 min.
Working temperature:	approx. 22°C (Room temperature)
Initial setting time:	9 - 11 min.
Setting expansion:	approx. 2.20% (type 100 liquid, 100%)

Conventional preheating:

If a conventional preheating is necessary should allow the mould to set for at least 60 minutes, then place it in a cold furnace and heat to the required temperature at a heat rate of 5 - 7°C per minute. Holding stages are not required when using Presto Vest II.

Concentration of the expansion liquid:

24 ml liquid / 100 g powder

High-gold-content precious metals:

Liquid ratio type 100 / distilled water: 55 % / 45 %

100 g	160 g
13 ml / 11 ml	21 ml / 17 ml

Semi-precious metals:

Liquid ratio type 100 / distilled water: 60 % / 40 %

100 g	160 g
14 ml / 10 ml	23 ml / 15 ml

Bonding alloys:

Liquid ratio type 100 / distilled water: 65 % / 35 %

100 g	160 g
16 ml / 8 ml	25 ml / 13 ml

Semi-precious bonding alloys:

Liquid ratio type 100 / distilled water: 70 % / 30 %

100 g	160 g
17 ml / 7 ml	27 ml / 11 ml

Precious free alloys (CrCo):

Liquid ratio type 100 / distilled water: 85 % / 15 %

100 g	160 g
20 ml / 4 ml	32 ml / 6 ml

Devesting:

Single crowns and inlays can be devested just 30 minutes after casting, but bridges have to be left to cool to room temperature before devesting. Presto Vest II is easy to devest using plaster nippers and, if necessary, a pneumatic chisel. Soak moulds for a few minutes before devesting to minimise the amount of dust produced. Sandblast with aluminium oxide and/or glass beads at 3 bars.

Please note:

Technical recommendations are based on tests and findings from work in our development laboratory and can only be regarded as guidelines. SILADENT products are subjected to strict quality controls. We reserve the right to make technical changes.

Warning:

Investments contain quartz and cristobalite. Do not breathe in dust. Danger of lung diseases (silicosis or cancer). Use a suitable dust mask!

Further questions: Contact SILADENT technical services (Tel.: +49 (0) 53 21-37 79 25/26) or our sales representatives.

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