# Silaform<sup>®</sup> Gingiva 1:1

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# Silaform<sup>®</sup> Gingiva 1:1

## Instructions for use

**Application:** Silaform Gingiva soft (REF 127300) is an addition curing silicone for the production of gingival reproduction.

<b>Technical data:</b>	Mixing ratio A : B:	1 : 1 Automix
	Working time at 23 °C:	1 min. 30 sec.
	Curing time at 23 °C:	7 - 8 min.
	Shore A hardness:	70

**Preparing the cartridge:** Place the cartridge in a 50 ml dispenser (fits all commercially available dispensers). Remove and dispose of the cartridge seal. The seal must not be reused, as the cap may contaminate the material if it is refitted incorrectly. When using initially and before fitting the mixing tip, extrude a small amount of material until a uniform consistency of silicone is discharged from the two outlets. Fit the mixing tip onto the cartridge, rotate to lock it in position and then dispense as required. After use, the mixing tip should remain fitted on the cartridge as a seal until the next application.

**Direct fabrication technique, application in the impression:** Spray Silaform Gingiva Sep (REF 127301) onto the area of the impression to be reproduced to insulate the area. Then apply Silaform Gingiva into the impression using uniform pressure. The mixing tip should remain in the material during application to avoid air bubbles. After curing, adapt the gingival mask using a scalpel if required and fabricate the stone model in the usual way.

**Indirect fabrication technique on the model:** First, fabricate the stone model in the usual way. Fabricate a silicone overcast using an addition-curing impression material (e.g. Silaform 90 extra hard, REF 102704), which should cover all the areas of the model that are to be reproduced. Alternatively, an existing impression can also be used as an overcast following removal from the impression tray. After removal, cut away any areas in the overcast or impression outside the gingival mask region that interfere with the fit to ensure easy repositioning on the model. Prepare the areas of the model that are to be reproduced with Silaform Gingiva soft. Insert one or more syringe inlet channels into the silicone cast (palatally or lingually). Additional ventilation channels should be inserted opposite to the syringe inlets. Clean the overcast and apply Silaform Gingiva Sep (REF 127301) to the inside to ensure separation. Replace the overcast on the model, check that it fits exactly and then fix it in position. Place the mixing tip into the inlet opening and syringe the material into the mould using uniform pressure until material exudes from the ventilation channels. After the gingival mask has cured, carefully lift the silicone overcast, separate the material in the channels using a scalpel and then remove the overcast completely. The gingival mask can now be removed from the model and prepared using scissors, scalpel or rotary instruments.

**Important information:** Do not store above 25 °C. Do not use after the expiry date. Do not use latex gloves, as they can impair the setting reaction of the polyvinyl due to contamination.

**Toxicology:** Non-toxic, no absorption, but may cause irritation to the eyes. The product can be used without any risks to health, provided it is handled correctly and in accordance with standard occupational hygiene measures.

**Ecological information:** Current findings indicate that no adverse effects on the environment are to be anticipated, provided the product is used correctly and in accordance with current regulations. The product is non-biodegradable and insoluble in water when cured. It can be easily separated from water by filtering.

**Please note:** Cured silicone is not harmful to the environment and can be disposed of with domestic waste in accordance with local regulations. Our recommendations on technique are based on tests and findings from our development laboratory and should only be regarded as guidelines. SILADENT products are subjected to strict quality controls. We reserve the right to make technical changes.

Date of Information : 03/2012