Silavest Evolution

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DIN EN ISO 15912

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Working instruction

DIN EN ISO 15912

Description: **Silavest Evolution** is a phosphate-bonded, graphite-free precious investment speed casting investment for crowns and bridges in the speed casting technique with extraordinary working properties. **Silavest Evolution** is especially designed for casting BÄR-LIGHT® (selection GmbH -dental-) with a reduced preheating temperature of only 700° C with a reduced residence time in the preheating furnace (approx. 30 minutes). **Silavest Evolution** is also usable for all other non precious alloys with a preheating time of max. 850° C.

Technical data:	Mixing liquid	Expansion liquid Typ 100
	Mixing ratio (powder : liquid)	100 g– 22 ml Liquid
	Mixing under vacuum	60 seconds
	Working time	5 – 7 min.
	Working temperature	17°C - 19°C
	Initial setting time	11 – 13 min.
	Setting expansion	> 2.00 %

Waxing up: Optimum fit and surface quality are best achieved using patterns fabricated entirely from wax. Using resins often results in much poorer surface quality and cracks in the mould with the speed cast technique. This applies particularly to waxing up removable sections of attachment and telescopic crown work.

Crown and bridge work: Mix Non precious investment according to the instructions and pour it up to the patterns without vibration. Then vibrate gently into the patterns using a small instrument, if necessary. Fill the ring completely without any further vibration. Adhere strictly to the time given in the instructions for the speed technique.

Use of pressure: Investing under pressure is not necessary due to Non precious investment extremely homogenous grain distribution. If investing under pressure is preferred, the pressure should be released after 15 minutes so that it does not impede the setting expansion.

Attaching sprues: For crowns and bridges: the best method for casting channels is the double curve system or the one with the metal wire ring for wax with a diameter of 3.5 mm at most. The diameter of joints with the cast object must be Ø 2.0 mm. The distance between the cast object and the double curve/ring must be 3.0 - 3.5 mm depending on the diameter of the wire used. The cast objects must be arranged in the muffle so they are never in the center of heat. The thinner parts (e.g. the edges of the crowns) must point towards the wall of the muffle and must be at least 5.0 mm away from its edge, while the thicker parts must be a much greater distance away. Use muffle sizes 3 and 6 for preference.

Mixing ratio per 160 g sachet: 35 ml liquid

Liquid concentration (expansion liquid type 100)*:	Single crowns and smaller bridges:	65 %
	Secondary cone crowns:	70 % - 85 %
	Telescopic secondary crowns:	90 % - 100 %

* The different liquid concentrations depends to the different contractions of the used modeling wax and other modeling materials. Depending to the situation you must adjust the liquid concentration.

Speed casting: **Silavest Evolution** is preferably preheated in the speed casting method. Preheat the furnace to 700 °C and, 20 minutes after mixing the investment, place the mould in the furnace with the sprue hole facing downwards. Adhere strictly to the prescribed time, which is measured from the first contact of the powder and liquid. Moulds should not be placed directly onto the floor of the furnace. We recommend using either a grooved tray for collecting wax or special stilts from a pottery craftshop. Caution! During the first 15 minutes the furnace fan and extractor should remain switched off and the furnace door should not be opened because of the risk of instant combustion. If a preheating temperature other than 850 °C is required, the temperature can be adjusted 15 minutes after placing the mould in the furnace. Moulds are ready to cast after heat soaking for 30 - 40 min. depending on the size of the ring.

Recommendations: These technical hints are based on practical experience and experiments conducted in our R&D laboratory. The values given are purely indicative. All products are subjected to strict quality controls. We reserve the right to make any technical modifications considered necessary.

Warning: Investments contain quartz and cristobalite. Do not breathe in dust. Danger of lung diseases (silicosis or cancer). Use a suitable dust mask!

Date of the information: 02/2015

Printing date: 12.11.2015

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SILADENT Dr. Böhme & Schöps GmbH