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Splint PMMA BioStar

Instructions for Use

Composition: Polymethylmethacrylate (PMMA)

Technicle Data

Flexural strength	ISO 20795-2	MPa	64
Flexural modulus	ISO 20795-2	MPa	1950
Fracture toughness	ISO 20795-2	MPam ^{1/2}	2,0
Charpy impact strength, +23°C	ISO 179/1eU	kJ/m ²	45
Density	ISO 1183	g/cm ³	1,17
Water sorption	ISO 20795-2	µg/mm ³	< 24
Solubility	ISO 20795-2	µg/mm ³	< 0,13
Residual monomer content	ISO 20795-2	%	< 0,7

Indications

Splint PMMA BioStar are transparent dental milling blanks out of PMMA for the manufacture of splints, therapeutic splints and bite regulators for long-term application in the oral cavity up to 12 months.

Contraindications

- intolerance to the contained components
- crowns, bridges, model casting with clasps, bars, narrow sublingual bars and transversal connectors

Safety information

Please pay attention to the information in the material safety data sheet in its current version. Avoid the inhalation of dust particles during processing. Wear protection gloves, safety glasses and a dust mask to prevent irritations of eyes, skin and respiratory system.

Handling and storage

The medical device may only be processed by trained personnel. Only use clean, defect-free material from undamaged packaging and ensure that the blanks are stored at room temperature, in their original cartons, in a dry location and protected from light. Avoid contact with fluids.

Milling and further processing:

Milling must be carried out with a milling strategy adapted to PMMA. Only use one-edged milling tools in the milling machine or wet grinding procedures!

Further processing:

Tungsten carbide burs like the one you use for soft reline materials are suitable for rough adaption. For the finalization common cross-toothed millers can be used. Polishing should be carried out with regard to a low contact pressure, otherwise the resulting heat might lead to smearing or melting.

Cleaning

- ultrasonic bath (max. 40 °C) with mild cleaning agent
- Do not use the following methods for cleaning:
 - steam blasting
 - autoclave method (e.g. sterilization / disinfection)
 - organic solvents and acids (e.g. acetone, hydrochloric acid)
 - highly concentrated alcohols
 - chemical denture cleaners
 - abrasive cleaning agents

Disposal

Remains of PMMA blanks can be disposed of in the regular household garbage.

We reserve the right to make changes because of the continuous development of our products.

Please find the current version of the instruction manual on our website:
www.siladent.de

This version replaces all previous versions.

Date of information: 05/2025

Explanation of symbols

Manufacturer

Date of manufacture

Consult instructions for use

Use before

Batch code

Catalogue No.

Keep dry

Protect from direct sun light

Medical Device

Unique device identification

Distributor

Content Quantity

Caution: U.S. Federal law restricts this device to sale by or on the order of a dentist.

Printing date: 22.05.2025

Hersteller / Manufacturer:
Dental Direkt GmbH
Industriezentrum 106 – 108
32139 Spence, Germany

SILADENT Dr. Böhme & Schöps GmbH