V-Alloy II CE 0124

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V-Alloy II CE 0124

Instructions for use

Area of application: Non precious alloy based on cobalt for removable dentures acc. DIN EN ISO 22674, type 5.

| Composition: (in % by mass) | Technical data: (guidelines) |
|--|----------------------------------|
| Co 64,0 | Proof stress 0.2 % (MPa) 579 |
| Cr 29,0 | Elongation at rupture (%) 6,3 |
| Mo 6,0 | Modulus of elasticity (GPa) 211 |
| | Tensile strength (MPa) 725 |
| Other constituents: Si, C, Mn, Fe | Vickers hardness VH 10 386 |
| | Density (g/cm ³) 8.4 |
| | Melting range (°C) 1.350-1.406 |
| | Casting temperature (°C) 1.445 |

Tolerances in the composition of the alloy (% by mass) are within the range of current DIN standards.

Product description:

V-Alloy II is an alloy for denture bases is distinguished by an outstanding fluidity, which grants to fill even the thinnest details of the frame, Moreover its molecular structure allows to obtain smooth, compact surfaces with little oxide formation. It is highly corrosion resistant and is free of toxic elements like beryllium, indium and gallium ($\leq 0,1\%$).

Investing and casting:

V-Alloy II can be used with phosphate-bonded CrCo investments such as Granisit®, Micro, Jet 2000 or Granisit® XF Speed. Prepare the moulds for casting according to the SILADENT System. Preheating range is 850 °C – 1.000 °C. We recommend 920 °C. Heat soak at casting temperature for 30 minutes.

The alloy is cast with an induction machine when all the cylindrical ingots have melted together and a final shadow runs over the molten metal. Considerable experience is required when melting high-grade alloys automatically, and it is most important that the heating is set exactly. Accurate heat setting is essential to protect the properties of the alloy. Light the flame and adjust the pressure in order to have a blue nucleus about 4 mm long. Place some metal ingots into the preheated crucible. Move the flame over the crucible rotating the tip, in order to evenly heat the alloy. As soon as the ingots are wholly melted, and after disappearing of their shadows, start centrifugal action. Casting V-Alloy II repeatedly is not recommended. During re-use pay attention that the alloys are from the same batch. Mix old and new alloys in a ratio of 1:1.

Preparing and polishing:

Cool down the muffle to room temperature, Sandblast the framework with aluminum oxide 110-250 µm (3 - 4 bar) and continue the working cycle according to your personal method. Remove carefully any small casting blows, go over all the corners and edges with a round bur and then smooth the surfaces with rubber abrasives and polish them.

Soldering and welding:

Presolder using SILADENT CrCo solder (1mm, REF 102878, 2mm, REF 102807). Laser weld using a SILADENT CrCo laser (REF 102806) welding rod.

Safety hints:

Metal dust is hazardous to health. For finishing and sandblasting use a suitable extraction system and / or face mask. The dental alloy can affect MRI results. Each supply is identified by a lot-number. Take note of this number on the patient's file in order to trace down the product.

Contraindications, side effects and interactions with other dental alloys:

In cases of hypersensitivity (allergy) to the constituents of the alloy, discontinue its use. A Patch-Test is thus advisable. In individual cases, hypersensitivity reactions (allergies) and electrochemically induced local dysaesthesia have been reported, such as changes in taste and irritation of the oral mucosa. Galvanic effects can occur under proximal or antagonistic contact with dentures of different alloys.

Storage conditions:

No special storage conditions are required.

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