
Gebrauchsanweisung | Instructions for use | Mode d'emploi | Modo de empleo | Modalità d'uso

rematitan®

rematitan® powder 15 - 45 µm

CE 0483

DE

rematitan® powder 15 - 45 µm

Pulver

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EN

rematitan® powder 15 - 45 µm

TiAlV – dental alloy based on titanium in ELI quality (class 23) for firing and partial denture technique, Type 4 acc. to DIN EN ISO 22674 / DIN EN ISO 9693 / DIN EN ISO 5832-3 for processing with laser melting machines

Composition (% by mass)

Ti	Al	V
90,0	6,0	4,0

Additional elements < 1%: Fe, O

Tensile strength*	R _{p0,2}	950 MPa
Tensile strength*	R _m	1005 MPa
Hardness*	H	330 HV10
Elongation at rupture*	A ₅	10%
Modulus of elasticity*	E	130 GPa
Density	ρ	4.4 g / cm³
Solidus temperature	T _i	1605 °C / 2426 °F
Liquidus temperature	T _l	1650 °C / 2588 °F
CTE (25 °C - 500 °C / 77 °F - 932 °F)	α	10.1 x 10 ⁻⁶ K ⁻¹

* Variations of ± 10 % are possible depending on the batch.

The technical data represent typical values that are achieved when using rematitan® powder. A laser melting machine of the type Concept Laser M2 Series 5 was used for the production of test specimens. The above-mentioned values were achieved with a layer thickness of 30 µm and 60 µm.

Dear customer,

Thank you for choosing a quality product from Dentaurum. It is essential to read these instructions carefully and adhere to them to ensure safe, efficient use and ensure that you and your patients gain full benefit. Instructions for use cannot describe every eventuality and possible application. In case of questions or suggestions, please contact your local representative.

As our products are regularly upgraded, we recommend that you always carefully read the current Instructions for use supplied with the product and stored in the internet at www.dentaurum.com, even though you may frequently use the same product.

1. Manufacturer

Dentaurum GmbH & Co. KG | Turnstraße 31 | 75228 Ispringen | Germany

2. General product description

Non-precious alloy for the manufacture of medical products by means of an additive powder bed fusion process, e.g. selective laser melting.

3. Indication

rematitan® powders for additive manufacturing are dental materials based on titanium. rematitan® powder is particularly suitable for the manufacture of fixed and removable prosthetic restorations, appliances and metal-ceramic frameworks.

4. Digital modeling

The dental products are modeled using a CAD software which has been approved for this purpose. Construction must be in line with dental regulations in order to fulfil mechanical requirements. The finished product must, however, have a material thickness of at least 0.3 mm.

5. Processing in the powder bed system

Before processing in the powder bed system begins, it must be clarified that the system is suitable for processing the material. Please contact your machine manufacturer in order to clarify this or contact us. We will be pleased to help you check the system requirements.

6. Post treatment

Powder residues on the components manufactured by the additive manufacturing process must be removed with due care, taking national safety precautions into consideration.

7. Stress relief heat treatment

Stress relief heat treatment for products manufactured by additive manufacturing technology ensures that the material properties are as indicated and guarantees that the components are of high quality and lasting durability. If the heat treatment program is altered, results may differ. This, in turn, may void the guarantee.

Stress relief heat treatment optimizes the material structure of the titanium alloy. Inner stress is relieved and the fit of the dental component is improved. The heat treatment of the components is recommended with a build plate in a suitable furnace exclusively with exclusion of oxygen and thus under a shielding gas atmosphere or in an evacuated environment. When using a gas box, observe the gas flow rates recommended by the manufacturer. Argon is preferably used as the shielding gas. Heat treatment under ambient air must be avoided.

Heat treatment classic:

- Heat furnace including build plate with supported parts within 4 hours to 850 °C/1562 °F.
- Hold the temperature at 850 °C/1562 °F for 90 minutes, then switch off the heating. If a gas box is used, it can be removed from the furnace. Ensure flushing with shielding gas continues.
- Remove the build plate once room temperature has been reached.
- Processing of the components can then be continued.

8. Post treatment

The components can be removed from the build plate using a separating disc or band saw. Support residues can be removed and the surface smoothed using a finely-toothed hard metal bur. It is important to observe the warnings and precautions to be taken.

9. Powder processing

Excess powder can be passed through a sieve with a mesh of approx. 63 µm. The safety-relevant aspects must be observed in accordance with the applicable national regulations. The separated powder should undergo visual inspection and can, if considered in good order, be further processed in the powder bed system. Alternatively it can be stored, dry and under seal, in its container.

10. Welding

Suitable welding processes such as laser welding produce mechanically strong and corrosion-resistant joints. It is important to observe the framework design, the surface structures, the welding sequence and the welding parameters recommended for each individual laser machine. Suitable welding wires are available as filler material for all titanium materials, e.g. rematitan® wire on coils 0.4 mm (REF 528-039-50) or 0.7 mm (REF 528-040-50) and rematitan® straight wire 1.0 mm (REF 528-041-00) or 1.2 mm (REF 528-042-00).

11. Ceramic veneering

Titanium materials for the bonding technique can be veneered with all suitable ceramic materials, e.g. ceraMotion® Ti (Dentaurum), as long as the CTE is suitable for the material. The highest recommended firing temperature is 880 °C / 1616 °F. The ceramic should cool down over a long period in the furnace unless the ceramic manufacturer has other specifications.

12. Contraindications and adverse reactions

Signs of intolerance to titanium and titanium alloys are extremely rare if the manufacturer's Instructions for use are adhered to. If the patient has a proven allergy against any component within the alloy, this alloy must not be used for safety reasons. There have been individual reportings of local irritations which were electrochemically induced. If various alloys have been used, it is possible that galvanic effects may occur.

Exposure to dust and vapors may cause irritation of the eyes and/or respiratory tract. Please adhere to warnings and safety data sheets.

13. Warnings and precautions

Metal shavings and metal dust may be produced as the alloy is mechanically processed. The mechanical process will also create heat and possibly burs. Processed material may therefore be hot and/or have sharp edges. It is possible that dust may be generated when the package is opened or the rematitan® powder is transferred. Metal dust is combustible. Any type of ignition source should be kept at a distance. Due to the oxygen affinity of titanium, there is a risk of burns and deflagration, particularly when the material is hot. Metal fire powder extinguishers should always be positioned in close proximity to extinguish a fire if necessary. They should be in working order. The workpiece must be cooled during mechanical processing in the presence of oxygen. Processing at temperatures above the solidus temperature may generate vapors. The workplace should be fitted with a suitable suction or ventilation system. We generally recommend that gloves, protective clothing, safety glasses and respiratory protection equipment are worn. Should there be contact with the skin, please wash the affected area with soap and water for a few minutes. Should the product come into contact with the eye, please rinse the eye with a gentle stream of water or saline solution for at least 15 minutes.

Please observe the relevant safety data sheet (see www.dentaurum.com).

Interproximal or occlusal contact with other metallic elements and the use of different types of alloys in the same oral cavity should be avoided to prevent galvanic reactions or local irritations in individual cases. The alloy should not be ground or polished in the patient's mouth. There is no scientific evidence on the safety or efficacy of treatment of pregnant women or nursing mothers or children.

14. Reuse

The material may be changed in its composition and properties if reused. In principle, we recommend using new material to manufacture high-quality prosthetic restorations.

15. Conditions for storing

Normal changes to the surrounding conditions (e.g. temperature, pressure, light) have no negative impact on the product characteristics of titanium and titanium alloys. rematitan® powder should be stored in a dry place in its original packaging with the lid closed tightly.

16. Disposal

Please adhere strictly to the rules that apply for your region, bearing in mind the details outlined in the safety data sheets.

17. General information

Should the user and/or the patient become aware of serious problems arising from the use of the product, it is important that the manufacturer and the competent authority in the country in which the user and/or the patient is resident is informed accordingly.

For more information on our products and services, please visit www.dentaurum.com.

In case of questions or suggestions, please contact your local representative.

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