## CoCr BioStar

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CE 0124

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**Area of application:** CoCr BioStar is a beryllium- and nickel-free non-precious bonding alloy for the production of milled crowns and bridges. Only for professional user!

#### Composition (in % by mass):

Co: 61.65 Cr: 27.75 W: 8.45 Si: 1.61

**Other constituents:** Si; Mn

Tolerances in the composition of the alloy (% by mass) are within the range of current DIN standards.

#### Technical data (guidelines):

Proof stress (MPa):	440
Elongation at rupture (%):	14.7
Modulus of elasticity (GPa):	240
Coefficient of expansion:	25-500 °C $14.5 \times 10^{-6} K^{-1}$
Tensile strength (MPa):	819
Vickers hardness VH 10/30:	275
Density (g/cm <sup>3</sup> ):	8.55
Melting range (°C):	1.395 -1.418
Max. oxide firing temperature (°C)	980

**General guidelines for handling:** This instruction for use includes important processing steps and recommendations for CoCr BioStar.

**Indication:** Crowns and wide-span bridges up to 16 units (max. 4 pontics in the anterior region and max. 3 pontics in the posterior region), implant retained superstructures and implant bars.

**Milling:** CoCr BioStar is construed for CNC milling machines. Please follow the instructions and parameters of the respective manufacturer of CAM Software and the CNC milling machine.

**CAD:** The design should be done with appropriate CAD software. Please consider an anatomically reduced framework design for the veneering with ceramic. The wall thickness should not be less than 0.3 mm. Choose a sufficient connector dimension (6 - 9 mm<sup>2</sup>). Sharp edges and undercuts should be avoided.

**Cutting out the frameworks from the blank:** Remove the milled frameworks with suitable cutting tools and smoothing the supports.

**Soldering / Laser welding:** CoCr BioStar can be soldered with all suitable solder. CoCr BioStar parts should not be soldered with gold or palladium solder. CoCr BioStar is also ideally suitable for laser-welding.

**Preparation before ceramic veneering:** The frameworks can be elaborated with standard carbide cutters, look for smooth transitions and avoid overlapping material. Please use the same cutter for one alloy to avoid contamination. The minimum thickness of the prepared coping should not be less than 0.3 mm. It's recommended to sandblast the frames with minimum 110 µm of Aluminium oxide with 3-4 bar and clean with steam cleaner. Oxide firing is not mandatory but can be done as an option for 5 minutes at 980°C (cleaning firing). The frame needs to be sandblasted again to remove the present oxide layer thoroughly. In the end the cleaning by steam cleaner is mandatory. If you use a ceramic bonder please consider the instruction for use of the manufacturer.

**Handling conditions / Safety:** Metal dust is harmful to health. Use when grinding and sandblasting dust extraction and respirator with filter FFP3 – DIN EN 149.

**Contraindications and side effects:** If the instructions are observed during the production processes, incompatibilities with non-precious dental alloys are extremely rare. In case of a proven allergy against an ingredient of this alloy, the alloy must not be used for safety reasons. In exceptional cases, electrochemically induced, local irritations have been reported. When different alloy groups are used, galvanic effects might occur. Please inform your dentist regarding the contra-indications and side effects. Any serious incident that involve the product must be reported to the manufacturer and the competent authority in the accorded country.

**Single-use:** Used milling discs should be not processed for further use as a medical device.

**Disposal Instructions:** Consult the material safety data sheets or national regulations for disposal. Dispose of CoCr BioStar residues and dust in an environmentally friendly manner. Grinding dust must not enter groundwater, water bodies or sewers. Address waste exchanges for recycling.

**Storage conditions:** Temperature, humidity or light has no effect on the product properties.

**Further questions:** Contact our technical services (Tel.: +49 (0) 53 21 – 37 79 25/26) or our sales representatives.

Date of information: 12/2018

Printing Date: 13.12.2018