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Stand der Information:  
09/2020

## Keralloy® N Instructions for use

**Area of application:** Non precious casting alloy based on nickel for metal ceramic accord. DIN EN ISO 22674, type 3.

Composition:	(in % by mass)	Technicle data:	(guidelines)
Ni	63,0	Proof stress (MPa)	360
Cr	25,0	Elongation at rupture (%)	26
Mo	9,0	Modulus of elasticity (GPa)	191
Si	2,0	Coefficient of expansion:	
Nb	1,0	25 - 500 °C	13,7 x 10 <sup>-6</sup> K <sup>-1</sup>
		25 - 600 °C	14,0 x 10 <sup>-6</sup> K <sup>-1</sup>
		Tensile strength (MPa)	602
		Vickers hardness HV 10	180
		Density (g/cm <sup>3</sup> )	8,3
		Melting range (°C)	1.298 -1.344
		Max. oxide firing temperature (°C)	950
		Casting temperature (°C)	1.400

Tolerances in the composition of the alloy (% by mass) are within the range of current DIN standards.

### Product description:

Keralloy® N is a non-precious, nickel based bonding alloy. Modeling with a minimum thickness of 0.3mm on simple structures and 0.5mm in complex cases or with patients with bruxism. Avoid elbows. The structure of the alloy generates less oxidation also after repeated casting and offers a gut processing and polishing. Its thermal expansion coefficient is ideal for every ceramic of the latest generation. Keralloy® N is highly corrosion resistant and does not contain any beryllium, indium or gallium ( $\leq 0,1\%$ ).

**⚠** This alloy contains nickel and should not be processed by individuals with known nickel sensitivity (allergy). Metal dust containing nickel should not be inhaled. Keralloy® N should not be used with anyone allergic to nickel.

### Investing and casting:

Keralloy® N can be used with phosphate-bonded crown and bridge investments, such as Premium and Presto Vest II. It is preheated to 850 - 950 °C and heat soaked for at least 30 minutes before casting. Follow the recommendations in the SILADENT technique instructions when attaching sprues. Use a separate ceramic crucible for casting Keralloy N high-grade, non-precious alloy to prevent it being contaminated by other alloys. To avert shrinkage blisters and cracks avoid hollows, sudden elbows and impurities. Wax walls should be 0.5 mm thick, so that after manufacturing process metal wall thickness is at least 0.3 mm. Transitions from wax to metal should be realized with gradual passages. High-frequency melting equipments: When last ingots have melted in the preheated ceramic crucible and incandescence shadow has disappeared, start melting process immediately. Do not use any flux. Open-flame melting: place ingots in the pre-heated ceramic crucible and heat them evenly with circular movements. When ingots have melted, start centrifugal unit. Use multi-flame welding torches only. Do not use any flux. Sample values for flame regulation: acetylene 0.4 bar / oxygen 2 bar ; propane 0.2 bar / oxygen 2 bar; leading pressure methane / oxygen 2 bar. Overheating the melted material could cause shrinking hollows, microporosity and building of a coarse-grained structure, and cause bridges to break or ceramic investments to crack. Casting Keralloy N repeatedly is not recommended, as this could alter both its composition and properties. It is nonetheless possible to re-use sprues and cones only once, given that new metal coming from the same lot is added in equal quantity.

### Manufacturing and polishing:

Cool down the muffle at room temperature, remove the investment and sandblast with aluminum oxide 100 to 150 µm at a pressure of 3-4 bar. Continue manufacturing process using fraises for hard metal. It is recommended an homogenization and degassing treatment of the casted piece at 1000°C for 10 minutes in atmosphere.

### Firing porcelain:

1. The surface has to be sandblasted with disposable equipment in aluminium oxide of 110 to 150 µm at a max. pressure of 3-4 bar and then steamclean. After cleaning, the framework should not be touched anymore.
2. Fire for 5 min. at 950°-980°C.
3. After firing, the oxide layer has to be carefully sandblasted using disposable aluminium oxide equipment 110 to 150 µm at a pressure of 2.5 - 4 bar, then steam clean.
4. If the framework has a homogeneous grey surface, then you can apply and fire ceramic, following its manufacturer's instructions.
5. It is possible to use all ceramics normally available on the market for NiCr alloys. Apply manufacturing instructions of the ceramic producer (e.g. VITA VM13, Shofu Vintage Halo). Slow cooling is recommended. It's possible to use Bonding.

### Soldering and welding:

Presolder using SILADENT CoNiCr solder (REF 102877), Laser weld using a SILADENT Cr-Co laser welding rod (REF 102806).

### Safety hints:

Metal dust is harmful to your health. When deflasking and blasting use a suction extraction system and breathing mask! The dental alloy can affect MRI results. Each supply is identified by a lot-number. Take note of this number on the patient's file in order to trace down the product. It is recommended that the patient be made aware of the possibility for dental alloys to affect MRI results.

### Contraindications, side effects and interactions with other dental alloys:

In cases of hypersensitivity (allergy) to the constituents of the alloy (a patch-test is thus advisable). Hypersensitivity reactions (allergies) or electrochemically-induced local dysaesthesia have been reported in individual cases. Galvanic effects can occur under proximal or antagonistic contact with dentures of different alloys. If lasting, electrochemically induced, local dysaesthesia arises from contact with other alloys, the restorations must be replaced with other materials.

### Storage conditions:

No special storage conditions are required.

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Mesa Italia Srl. packaging carries the mark

**CE 0123**

**Further questions:** Contact SILADENT technical services  
(Tel.: +49 (0) 53 21-37 79 25 / 26)  
or our sales representatives.

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**Commerciale:** SILADENT Dr. Böhme & Schöps GmbH, Im Klei 26, DE 38644 Goslar
**Fabbricante:** Mesa Italia Srl, Via dell' Artigianato, 35/37/39, IT 25039 Travagliato / Brescia
Mesa Italia Srl. segnare il confezionamento a