## TEK-1 Leg

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## TEK-1 Leg

### Instructions for Use

**Application: Non precious casting alloy based on cobalt for metal-ceramic accord. DIN EN ISO 22674, type 4.**

**Composition w%:**

Co 61.0% · Cr 27.0% · Mo 6.0% · W 5.0% · Si 1.0%

**Other constituents:** Mn, C, Fe

**Technicle data:**

Proof stress 0.2% (MPa)	475
Elongation at rupture (%)	8
Modulus of elasticity (GPa)	223
Coefficient of expansion 25°C-500°C	14.0 x 10 <sup>-6</sup> K <sup>-1</sup>
Vickers hardness VH 10	365
Density (g/cm <sup>3</sup> )	8.6
Melting range (°C)	1.350-1.385
Casting temperature (°C)	1.485
Max. oxide firing temperature (°C)	1.050

Tolerances in the composition of the alloy (% by mass) are within the range of current DIN standards.

**Investing and casting:**

TEK-1 LEG can be used with phosphate-bonded, graphite-free crown and bridge investments, such as Premium, Presto Vest II and TeleVest and for the TEK-1 one piece casting technique (TEK-1 Vest) as well. It is preheated to 850 °C – 920 °C, the holding time at end temperature - depends to the muffle size – until 45 minutes before casting. Follow the recommendations in the SILADENT technique instructions when attaching sprues. Use a separate ceramic crucible for casting TEK-1 LEG, to prevent it being contaminated by other alloys. Glaze the crucible once before first use. When using high frequency inductive melting, the casting process should be initiated as soon as the alloys surface tears up. When using a blow torch, crucible should be preheated. Place alloy ingots contacting each other, the better the melt is hooked up. Appropriate gas mixtures are propane/oxygen or acetylene/oxygen. Propane/compressed air are not appropriate, since it does not reach the necessary temperature (follow instructions on blow torch). Configure blow torch following engraved instructions on the torch (e.g. propane 0.5 bar and oxygen 2.0 bar). Only melt alloys with a reducing flame, all other configurations will harm the alloy. The flame is to be held in a distance of 7 cm from the alloy without rotation, while the cast remains permanently covered (otherwise inflow from the atmosphere). After the first shadows disappear from the cast, initiate casting before the oxide skin tears.

**Prepare and polishing:**

Cool down the muffle to hand temperature, do not knock against the casted model or against the casting cone during divesting. Risk of distortion! Sandblasting with pure aluminium oxide (250 µm) with 3.5 bar. Frames can be worked on with cross cut tungsten carbide burs or sintered stone. For polishing, of TEK-1 Leg frames use medium soft rubber wheels and cylinders of fine grain. After sandblasting, prepare the metal surface with corundum, diamond or carbide milling tools. For polishing of TEK-1 LEG don't use fine corns sized rubbers (wheels or roller). Afterwards polish the surfaces with a goat hair brush, for example, chewing surface should be polished with metal brushes under use of a suitable polishing paste. You will reach a gloss finish with linen or cotton buff under usage of high gloss polishing paste or similar polishing materials.

**Porcelain firing:**

Please follow the instructions of the porcelain manufacturer (e.g. Vita VM 13). Finishing glaze firings should be dealt with like preceding dentin firings, but without vacuum. The glaze firing should reach the desired surface gloss. It's possible to vary the end temperature up and down.

**Soldering and welding:**

A laser welder will have to be adjusted and configured individually. If the focus is variable, adjust it at between 1.1 and 1.3 mm. The pulse duration should be adjusted at between 11 and 13 ms, while tension should be adapted to pulse duration. Before welding, the chamber should be flooded with argon for at least 30 seconds. The ideal distance for the welding seam to be flooded with argon is 1 cm. Stained welding areas indicate an insufficient argon cover or a too high level of combined energy. Crack formation indicates too much energy or too long exposure time of the laser beam. While soldering, the parts must be prepared with a staggered tooth system cutter, without sandblasting. The soldering gap should not be more than 0,1 mm at the fixed assembly parts (risk of distortion). If required a piece of original alloy can be fit into place and spot-tacked. Construct a refractory model (eventually using ceramic pins) and place it after setting on a refractory block. After setting, heat refractory model slowly to dry the object carefully (please consider torch settings, use gas/oxygen). Heat soldering object up to flow temperature of the solder, apply flux and proceed soldering (solder forming spheres indicates low temperature).

**Contraindications, side effects and interactions with other dental alloys:**

Incompatibility reactions are very seldom, when working instructions are strictly considered. In cases of hypersensitivity (allergy) to the constituents of the alloy, discontinue its use. In individual cases, electrochemically induced local dysaesthesia has been reported. Galvanic effects can occur under proximal or antagonistic contact with dentures of different alloys.

**Storage conditions:**

No special storage conditions are required.

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**Manufacturer:** The Argen Corporation · 5855 Oberlin Drive, San Diego · California 92121, USA

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Printing date: 10/2020

**Further questions:**

Contact SILADENT technical services  
(Tel.: +49 (0) 53 21 – 37 79 19)  
or our sales representatives.

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**Distribuidor:** SILADENT Dr. Böhme & Schöps GmbH, Im Klei 26, DE 38644 Goslar
**Fabricante:** The Argen Corporation · 5855 Oberlin Drive, San Diego · California 92121, USA
The Argen Corporation rotula el envase comercial con **CE 2797**
**SILADENT Dr. Böhme & Schöps GmbH**