# TEK-1 Sil

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**Typ:**  
Additionsvernetzendes  
1 : 1 Dubliersilikon  
DIN EN ISO 14356,  
Typ 2 - irreversible  
Dubliermasse

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### Technical data:

Mixing ratio A (white): B (white):	1 : 1
Mixing under vacuum:	40 seconds
Working time at 23°C:	> 5 min.
Curing time at 23°C:	30 - 45 min.
Tensile strength:	approx. 3.0 MPa
Elongation at break:	approx. 360%
Tear strength:	approx. 5.0 N/mm
Shore A hardness after 60 minutes:	12 - 15

### Application:

TEK-1 Sil is suitable for all duplicating flasking systems. Models to be duplicated should be at room temperature and should not be soaked. Any dirt on the surface is removed by spraying with Neutralit and by blow-drying carefully. Materials, such as gypsum hardener, grease, oils and some separators, can produce a negative reaction with TEK-1 Sil (refer to „Incompatible materials“ below).

When selecting the flask size, ensure there is a gap of at least 10 mm between the model and the flask so that the TEK-1 Sil can flow round the model. Mix TEK-1 Sil under vacuum and pour in a thin flow until the model is completely covered. The stabilising or base plate is then placed in position. TEK-1 Sil can be measured using the SILADENT 1 : 1 dispenser DA 2000 (REF 111503) because both components have almost the same viscosity.

### Chemical characterisation:

Polymerised dimethyl polysiloxane filler mixture.

### Dangerous reactions:

No known dangerous reactions, if stored and used properly.

### Toxicity:

None known

### Directions:

Store containers tightly sealed in a well-ventilated place. Keep containers dry. Keep away from water, alkalis, strong acids and oxidising agents.

### Ecological information:

Present experience suggests that there is no detrimental effect on the environment, if it is used properly and current regulations are observed. It is non-biodegradable and insoluble in water when cured. It can be easily separated from water by filtering.

### Please note:

Cured silicone is not harmful to the environment and can be disposed of along with domestic waste, as long as this does not contravene local by-laws.

Technical recommendations are based on tests and findings from work in our development laboratory and can only be regarded as guidelines. SILADENT products are subjected to strict quality controls. We reserve the right to make technical changes.

### Incompatible materials:

The following materials can produce a negative reaction with addition-curing silicones:

soft polyvinyl chloride, polysulphide MIL-S-8516, Mystik® 6207 adhesive tape, Mystik® 6215 adhesive tape, Scotch® cellophane strips, Scotch® 360 adhesive tape, Permacel® masking tape, Pliobond® adhesive, Neoprene rubber, Buna N rubber, Vitron® A rubber, colophony-based soldering flux, acid-based soldering flux, disinfectant, ISOLIT, cyanoacrylate glue, amino-hardened epoxy, Humiseal® 6215 adhesive tape, vinyl insulating tape, latex vacuum hoses, GRS rubber, natural gum, certain types of silicone and RTV rubber, SUPER-SEP, milling oils, sulphur compounds such as thiols, sulphides, sulphates, sulphites, thiourea and nitrogen compounds such as amines, amides, imides, acids.

If there is any doubt, we recommend testing the material before use.

**Further questions:** Contact SILADENT technical services  
(Tel.: +49 (0) 53 21 - 37 79 25 / 26)  
or our sales representatives.

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