TeleVest

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DIN EN ISO 15912, Typ 1, Klasse 2

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DIN EN ISO 15912, type 1, class 2

Description:

TeleVest is a low dust, phosphate-bonded, graphite-free precision speed cast investment material. It is suitable for casting all types of dental alloy and can be preheated up to 1,100°C. It is for the use when casting all types of telescopic crown (e.g. TeleRing-technique) and crown- & bridge technique with precious and non-precious alloys.

Waxing up:

The highest precision of fit and best surfaces are achieved by waxing up the pattern with wax alone or other standard pattern resins. For the TeleRing-technique the use of standard pattern resins is recommended.

Application:

Mix TeleVest according to the instructions and pour it up to the patterns without vibration. Then vibrate gently into the patterns using a small instrument, if necessary. Fill the ring completely without any further vibration. Adhere strictly to the time given in the instructions for the speed technique.

Use of pressure:

Investing under pressure is not necessary due to TeleVest's extremely homogenous grain distribution. If investing under pressure is preferred, the pressure should be released after 15 minutes so that it does not impede the setting expansion.

Attaching sprues:

When using the TeleRing-technique, it is preferable to attach 2 - 3 mm Ø wax sprues directly. When casting outer full copings, it is advisable to use Perawax reservoir sprues (REF 103203-05). Position the patterns in the casting ring, ensuring that they are not in the thermal zone. The thin sections of the patterns should point toward the side or top of the ring. The patterns should be 5,0 mm from the side or top of the ring - size 6 casting rings are preferred. If plastic sprues are used despite their known disadvantages, they must be coated with wax to prevent the mould cracking.

Concentration of the expansion liquid for C&B-technique:

High-gold-content precious metals	55 - 60%
Semi-precious metals	60 - 65%
Bonding alloys	65 - 75%
Semi-precious bonding alloys	70 - 80%
Non-precious alloys	90 - 100%

Technical data:

Mixing liquid:	SILADENT type 100 expansion liquid
Mixing ratio Powder/Liquid:	C&B-technique 100 g : 24 ml secondary crowns 100 g : 22 ml
Mixing under vacuum:	60 sec.
Working time:	5 - 7 min.
Working temperature:	approx. 22 °C (Room temperature)
Initial setting time:	9 - 11 min.
Setting expansion:	approx. > 2.20% (100% undiluted type 100 liquid)

Expansion mixing liquid - Concentrations TeleRing- technique using light-curing wax or pattern resin:

	0°	2°	4°	6°
Non-precious alloys	100%	100%	95%	90%
Precious alloy, high gold content	50%	50%	50%	45%
Semi-precious alloy	60%	60%	50%	45%
Metal-ceramic alloy, high gold content	50%	50%	50%	45%
Semi-precious metal-ceramic alloy	60%	55%	55%	55%

Concentration of the expansion liquid for secondary structures outer copings made with light-curing wax, Adapter or pattern resin:

	0°	2°	4°	6°
Non-precious alloys	-	100%	100%	90%
Precious alloy, high gold content	60%	55%	50%	45%
Semi-precious alloy	60%	55%	50%	50%
Metal-ceramic alloy, high gold content	60%	50%	50%	50%
Semi-precious metal-ceramic alloy	65%	60%	60%	60%

Please note:

For wax patterns, the concentration of the expansion mixing liquid should be reduced by 5 - 10%.

Speed casting:

Preheat the furnace to 850°C and, exactly 15 minutes after mixing the investment, place the mould in the furnace with the sprue hole facing downwards. Adhere strictly to the prescribed time, which is measured from the first contact of the powder and liquid. Moulds should not be placed directly onto the floor of the furnace. We recommend using either a grooved tray for collecting wax or special stilts from a pottery craftshop.

Caution! During the first 15 minutes the furnace fan and extractor should remain switched off and the furnace door should not be opened because of the risk of instant combustion. If a preheating temperature other than 850°C is required, the temperature can be adjusted 15 minutes after placing the mould in the furnace. Moulds are ready to cast after heat soaking for 45 - 60 min. depending on the size of the ring.

Conventional preheating:

Allow the mould to set for at least 60 minutes, then place it in a cold furnace and heat to the required temperature at a heat rate of 5 - 7°C per minute. Holding stages are not required when using TeleVest.

Please note:

Technical recommendations are based on tests and findings from work in our development laboratory and can only be regarded as guidelines. SILADENT products are subjected to strict quality controls. We reserve the right to make technical changes.

Warning:

Investments contain quartz and cristobalite. Do not breathe in dust. Danger of lung diseases (silicosis or cancer). Use a suitable dust mask!

Further questions: Contact SILADENT technical services (Tel.: +49 (0) 53 21-37 79 25/26) or our sales representatives.

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