Know-How Brings Quality Back to Life

Instructions for use

Device: **YuDent™** Dental PEEK Block

3.....	English
6.....	Deutsch
9.....	Español
12.....	Français
15.....	Italiano
18.....	Português

Intended use

For the manufacture of frameworks used in the production of:

- i) Full and partial removable dentures and overdentures;
- ii) Anterior and posterior crowns, and posterior 3-unit bridges (maximum 1 pontic);
- iii) Implant supported suprastructures.

Contraindications

Not to be used in bridges with 2 or more pontics.

Precautions

- The minimum thickness of secondary crowns and matrix housings should not be less than 0.7 mm.
- The minimum thickness for “conventional” retention grids should not be less than 2 mm.

- Minimum size of horizontal retainer in lower jaw 2mm x 8 mm.
- Minimum size of horizontal band in upper jaw 2mm x 8 mm.

Material Information

YuDent™ Dental PEEK Blocks are manufactured from PEEK:

Form:	Solid
Color:	Grey
Odor:	Odorless
pH (Value):	Not known
Boiling Point (°C):	Not known
Melting Point (°C):	343
Flash Point (°C):	Not known
Auto Ignition Temperatur (°C):	595
Oxidising Properties:	Not applicable
Vapour Pressure (Pascal):	Not known
Density (g/ml):	~1.3
Solubility (Water):	Insoluble

Production

Dental PEEK Block can only be processed on suitable milling machines with appropriate milling programs. The milling head must be suitable for machining PEEK polymer. Yunyi (Beijing) Medical Device Co., Ltd. recommends suitable adaptation of the CAM software by the software manufacturer to allow milling of Dental PEEK Block. For an exact fit, a 5-axis milling machine of stable construction is recommended. Yunyi (Beijing) Medical Device Co., Ltd. also recommends cooling the milling head for more accurate milling and to prolong its useful life.

Prosthesis removal from the Dental PEEK Block

Remove the prosthesis carefully from the blank using a hard metal cutter, and pay particular attention to the edges of the prosthesis.

Finishing

For finishing work, hard metal cross-cut burs are recommended. Plain cross-cut burs or special fine-toothed burs are particularly suitable. Further fine finishing of horizontal bands and retainers can be carried out with a silicon polisher. Before high-gloss polishing, Yunyi (Beijing) Medical Device Co., Ltd. recommends brushing on the polishing motor.

Safety Information

While machining the Dental PEEK Block the following safety precautions are recommended:

- Dust mask or dust extraction

Storage information

The Dental PEEK Block should be stored in dry storage and avoid exposure to direct sunlight. The PEEK from which the Dental PEEK Block is made is stable and can be stored for

an extended period of time (10 year shelf life). It has a working temperature range from cryogenic up to 250°C and hence the storage temperature range for the Dental PEEK Block is any ambient temperature and humidity.

Additional Information

For additional information contact info@yunyimedical.com

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Yunyi (Beijing) Medical Device Co., LTD.
Shanghai International Holding Corp. GmbH (Europe)