## Zirkon BioStar

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CE 0124

gültig für Zirkon BioStar (Colour / S / L) / Zirkon BioStar Z / Zirkon BioStar S-HT/ HT Colour / HT ZZ

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### General description

Machinable blank made of yttrium stabilized, presintered zirconium dioxide for milled production of crowns and bridge frameworks on commercial CAD/CAM systems or hand-operated copy-milling machines. Ceramic material according to EN ISO 6872 Typ II, class 5.

### Indication

The different Zirkon BioStar material types are indicated for crowns and bridges (like caps, interlink and all-ceramic restoration with max 2 pontics at anterior and max 4 pontics at posterior bridges), up to bar constructions, telescopes, Maryland bridges, Inlays, Onlays and individual implant abutments.

### Contraindications

In the case of bruxism, inadequate hard tooth structure or insufficient space, the use of full ceramic crowns and bridges is not recommended or only advised in some cases. The same applies in the case of inadequate oral hygiene or inadequate coordination of masticatory movements. A full ceramic crown or bridge may also not be used in the event of unsatisfactory preparation results or in the case of production faults.

### Preparation information

The minimum wall strength values in the following table are merely reference values. The decision on the manner and undertaking of the preparation is the responsibility of the dental surgeon and dental technician (this includes the validation of the gluing process).

#### The wall strength should be at least:

Posterior teeth:	occlusal: 0.9 mm	circular 0.5 mm
Anterior teeth:	incisal/occlusal: 0.7 mm	circular 0.5 mm
Abutment crowns of bridge frameworks:	incisal/occlusal 1.0 mm	circular 0.7 mm

#### The connector surface should be at least:

Anterior bridge with one Pontic:	7 mm <sup>2</sup>
Anterior bridge with two pontics:	9 mm <sup>2</sup>
Posterior bridge with one Pontic:	9 mm <sup>2</sup>
Posterior bridge with two pontics:	12 mm <sup>2</sup>

### Milling

Milling can be undertaken with most commercially available CAD/CAM milling machines or hand-operated copy-milling machines. The machine manufacturer's manual should be consulted regarding the selection of the milling cutters and the necessary revolutions/minute. In order to avoid staining or a loss of strength, the blank should if possible not come into contact with liquids (water, oils) or fats (handcream) during milling and related handling of the blank. Wet machining of the dental blanks in their presintered state is in principle possible, but does require that the restoration be dried completely before sintering. Drying can be undertaken in a drying chamber at approx. 80° C or under an infrared lamp. Alternatively, a low-temperature firing program can be selected which will not damage the blank. Notches, sharp edges or slight damage, also on the underside of the framework, may reduce the strength of the completed sintered work. Consequently, it is already important to ensure that surfaces are well rounded and smooth when milling.

Zirkon BioStar HT Colour shows a higher hardness than the other qualities. Use a less feed motion, higher milling speed and synchronous milling is recommended.

**Inhalation of dusts may irritate the respiratory tract. Contact with the skin and eyes cause mechanical irritation.**

### Coloring

If required, the frameworks may be adjusted to the desired tooth color with the standard coloring solutions which are licensed for zirconium oxide ceramics. In order to achieve uniform results, it must be ensured that the framework is completely clean, greaseless and dry before coloring. The instructions of the manufacturers of the coloring solutions should be followed for the further work procedures.

### Sintering

Recommended sintering parameters (modifications of these recommended parameters, according to the framework, are possible at your own's authority).

	Zirkon BioStar (S / L / Colour) / Zirkon BioStar S (Colour) / Zirkon BioStar HT (S-HT / Colour / ZZ)	Zirkon BioStar Z
<b>Sintering temperature</b>	1.450 °C	1.500 °C
<b>Soaking time</b>	2,0 h	2,0 h
<b>Heating rate</b>	up to 600 K/h	up to 600 K/h
<b>Cooling rate</b>	up to 600 K/h	up to 600 K/h

The sintering of the preparations can be undertaken in all standard furnaces which are suitable for zirconium oxide ceramics. The manufacturer's operating instructions must be followed in all cases. The shrinking factor depends on the lot and it is noted on the label, the blank and the inspection certificate. We recommend sintering in a sintering crucible in which the restoration is embedded in suitable granular material. One suitable type of granules is ZrO<sub>2</sub> sintering spheres. When using sintering crucibles and granules, the framework is placed in the granules with the occlusal side downwards. The contact with other frameworks or sinter-bearers should be avoided. If sintering wires or pegs are used, it must be ensured that there is no contact with the cap wall area.

Dental frameworks and crowns made of zirconium oxide should be fired in a furnace which is only used for these products. If additionally restorations made of other ceramic materials are fired or if liquid ceramics are baked in the same furnace, this may impede the sintering behavior or result in localized discoloration. If the firing of different materials in the same furnace cannot be avoided, then an empty firing cycle at maximum temperature is recommended before firing zirconium oxide ceramics.

### Zirkon BioStar S - Additional information for the instructions for use

Instructions for the use of Zirkon BioStar S in your Sirona inLab<sup>®</sup> or inLab<sup>®</sup> MCXL system:

General information: Basically, you can use your Zirkon BioStar S blocks as other blocks for your Sirona<sup>®</sup> system. Programming of the Sirona inLab<sup>®</sup> or inLab<sup>®</sup> MCXL System:

Place the Zirkon BioStar S block as usual. Chose from the selection of blocks "VITA In-Ceram<sup>®</sup>YZ" the corresponding block size for the process. Enter the "Z-Code" indicated on the package completely, including "Z", as for instance: \*Z24912A. Inicie milling

VITA In-Ceram<sup>®</sup>YZ, Sirona inLab<sup>®</sup> and inLab<sup>®</sup> MCXL Systems are registered trademark of the manufacturers.

### Facing of the frameworks

Every commercial facing ceramic for frameworks made of zirconium oxide can be used. Please follow the instruction manual of the manufacturer.

### Finishing after final sintering

Any necessary finishing of the sintered constructions must only be done using water cooling during grinding operations, thus avoiding localized overheating and damage to the material.

### Storage information

If the dental blanks are stored in original packaging with regard to temperature and humidity ("cool and dry", namely 10°C to 35°C and 30 to 70% relative humidity), then no alteration of the material properties is to be expected, even if stored for longer periods. The blanks must not be stored in direct sunlight.

### Disposal of material residues

Remains of the blanks and, of course, also sintered parts made of Zirkon BioStar HT Smile materials are categorized as chemically completely inert and are not subject to any special disposal stipulations.

### Further questions

Please contact our technical services (Tel.: +49 (0) 53 21 – 3779-25) or our sales representatives.

## Technical data sheets for the different types of Zirkon BioStar

### General description

Machinable blank made of yttrium stabilized, presintered zirconium dioxide for milled production of crowns and bridge frameworks on commercial CAD/CAM systems or hand-operated copy-milling machines. Ceramic material according to EN ISO 6872 type II, class 5.

**Chemical Formula:**  $ZrO_2 / Y_2O_3 / Al_2O_3$

**Chemical Name:** Yttrium stabilized tetragonal zirconium dioxid

	Zirkon BioStar Z	Zirkon BioStar / Zirkon BioStar S / L	Zirkon BioStar Colour	Zirkon BioStar HT (S-HT/ZZ)	Zirkon BioStar HT Colour
<b>Colour</b>	white translucent	white opaque	opaque coloured	white high translucent	translucent coloured
<b>Chemical characteristics in wt %:</b>					
ZrO <sub>2</sub> / HfO <sub>2</sub> / Y <sub>2</sub> O <sub>3</sub>	> 99.0	> 99.0	> 99.0	> 99.0	> 99.0
Y <sub>2</sub> O <sub>3</sub>	5.15 ± 0.20	5.15 ± 0.20	5.15 ± 0.20	5.15 ± 0.20	5.35 ± 0.20
HfO <sub>2</sub>	< 5.0	< 5.0	< 5.0	< 5.0	< 5.0
Al <sub>2</sub> O <sub>3</sub>	< 0.1	0.15 – 0.35	0.15 – 0.35	< 0.1	< 0.1
Fe <sub>2</sub> O <sub>3</sub>	< 0.1	< 0.1	0.04 – 0.06	< 0.2	< 0.2
Na <sub>2</sub> O	< 0.04	< 0.04	< 0.04	< 0.04	< 0.04
Er <sub>2</sub> O <sub>2</sub>	-	-	-	< 0.1	< 0.1
Co <sub>3</sub> O <sub>4</sub>				< 0.1	< 0.1
<b>Physical material data after final sintering</b>					
Density (g/cm <sup>3</sup> )	> 6.05	> 6.05	> 6.00	> 6.00	> 6.00
Bending Strength (MPa)	> 800	> 1.000	> 800	> 800	> 800
Thermal expansion (10 <sup>-6</sup> 1/K)	10.4	10.4	10.4	10.4	10.4
Open porosity	keine	keine	keine	keine	keine

**Radio activity** Activity<sup>238</sup> U (Bq/g) < 0.2      DIN EN ISO 6872

### Biocompatibility and chemical solubility

Cytotoxicity ISO      inoffensive      DIN EN ISO 10993 / DIN EN ISO 9363  
 Chemical analysis      < 1 µg/cm<sup>2</sup>      DIN EN ISO 6872

The data indicated relates to test specimens from which it was obtained and is not unconditionally applicable to other forms of same material. The date must, therefore, be regarded as indicative only.

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