SilaPress® Vario

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SilaPress® Vario

Instructions for use

SilaPress® owes its popularity to its extensive processing range, which enables the user to work particularly comfortably.

Technical data:

Mixing ratio Powder/Liquid:	10 : 7
Swelling time approx.	30 s*
Casting phase approx.	3 min*
Plastic phase approx.	20 min*
Polymerization (2-3 bar, 45-50°C) approx.	15-20 min**

*The provided processing times refer to a material- and room temperature of 23°C [73,5°F]. Higher temperatures shorten, lower temperatures extend these times.
The polymerization time of 10-15 minutes refers to a previous working time of up to 20 minutes. The time in the pressure unit should last up to 35 minutes after mixing.

Indications:

- Completion of model cast dentures
- Total upper and lower dentures by using the pouring technique
- partial or total relinings (indirect method)
- Additions and repairs

Contraindications: Patient contact with unpolymerized material or any of its components is contraindicated. If the patient is allergic to one or more ingredients of the product, the product must not be used.

Material / Purpose: Denture base resin, autopolymer based on methyl-methacrylate in the form of powder and liquid. Classification according to ISO 20795-1 Type 2 Class 1 and according to Council Directive 93/42/EEC annex IX Class IIa for removable and for fixed dentures.

Product features:

- universal use
- medium processing time
- very good flow ability
- mucous membrane friendly due to a low amount of residual monomer
- absolutely color stable due to a catalyst system without tertiary amines
- Cadmium free
- evaluated and certified biocompatibility
- Phthalate and BPA free

SilaPress Vario does not adhere to high cross-linked synthetic teeth. Synthetic teeth should be prepared (grinded) and/or supplied with mechanical retentions.

Suggested mixing ratio: 10 g of powder with 7 g of liquid. Free dosing is possible, however conformity of the materials physical properties can only be guaranteed if the components are weighed.

Preparation: The arrangement of the prosthetic teeth in wax need to be fixed with a key made of plaster or silicone. Key, model and teeth are boiled out thoroughly to remove any wax residue. The plaster model needs to be well watered to assure that no air is pressed from the plaster into the acrylic during polymerization inside the pressure polymerization unit.

Isolation: Apply two thin layers of alginate based isolation to the model, with the second layer, it must be ensured that the first layer has not yet dried. Marmosep K always shows reliable results. Important: The alginate insulation must be well dried before first contact with the resin to prevent possible white discolouration on the base surfaces.

Processing: Place the required amount of liquid into the mixing cup, then add the according amount of powder. When dosing individually add just as much powder as needed to absorb the liquid. Stir thoroughly until the material reaches a homogenous consistency. After a swelling phase of approximately 30 seconds, the material is castable for 2 to 3 minutes. During this time the material is poured into the precasts. With beginning of the plastic phase, the material becomes steady, does not flow from the precasts and can be modeled for approximately 20 minutes.

Suggested polymerization: After a maximum of 20 minutes the polymerization is carried out under a pressure of 2-3 bar and a water temperature of 45-50°C [113-120°F] for about 15-20 minutes (up to 35 minutes after mixing).

After completion: Until incorporation the denture should be stored in water. It is recommended to water the denture for at least 24 hours to achieve an even higher resorption of the residual monomer (3,8% according to ISO 20795-1) and thus reduce irritations of taste.

Storage and safety measures: The material must not be used after expiration date. Storage temperature must not exceed 25°C [77°F] Avoid direct sunlight. Keep containers closed after use. Keep away from heat/sparks/open flames/hot surfaces. No smoking. Dispose of contents/ container in accordance with local recycling regulations. Highly flammable liquid and vapour. Causes skin irritation. May cause an allergic skin reaction. May cause respiratory irritation. Avoid breathing dust/fume/gas/mist/vapours/spray. Wear protective gloves/protective clothing/eye protection/face protection. If on skin: Wash with plenty of soap and water.

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