SilaPrint model2 385 nm

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Produktbeschreibung

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Instructions for use

Product description

SilaPrint model2 is a light curing resin for the production of dental models using 385 nm LED based stereolithography. For obligatory post curing of SilaPrint model2 built parts, the high performance light curing unit Otoflash G171 is recommended, which is optimal balanced to the production process.

Processing

- During the application of SilaPrint model2 we recommend wearing personal protective equipment f. i. suitable gloves, safety glasses etc.
- SilaPrint model2 should be well shaken approx. 1 hour before use.
- Pour SilaPrint model2 carefully into the pre-determined container of the production unit.
- Remove all bubbles with a cleaned object.
- Select the product specific parameter settings for your process.
- When building process is finished a direct post treatment is recommended.
- After the platform is taken up a dripping off time of approx. 10 minutes is recommended.
- Resin coated parts get clean with isopropanol (97 %) within a time of about 10 minutes. For a particular effective cleaning with isopropanol the use of explosion-protected ultrasonic units is recommended.
- Post-curing
Otoflash G171 with 2 x 3000 flashes under inert gas conditions (nitrogen), rotate components in between
- Impurity due to the construction or a break of the material and impurity due to operation mistakes cannot be excluded. Thanks to the low viscosity however it is possible to filtrate SilaPrint model2. It is recommended to regularly take out the container of the production unit, to homogenize and filtrate the content. Upcoming entrained bubbles get eliminated through a holding time of approx. 1 hour with SilaPrint model2.

To avoid damage of material quality does not expose the liquid material to light under any circumstances. Variations from the described manufacturing process may lead to modified mechanical characteristics and colour variations of the SilaPrint model2 material.

SILADENT Dr. Böhme & Schöps GmbH

Im Klei 26
DE-38644 Goslar

Tel.: +49 (0) 53 21 / 37 79 - 0
Fax: +49 (0) 53 21 / 38 96 32

info@siladent.de
www.siladent.de

SilaPrint model2 385 nm

Instructions for use

Possible dangers

Please note the safety instructions for processing SilaPrint model2:

Hazard statements: **H315** Causes skin irritation. **H319** Causes serious eye irritation. **H317** May cause an allergic skin reaction. **H413** May cause long lasting harmful effects to aquatic life.

Precautionary statements: **P261** Avoid breathing dust/fume/gas/mist/vapours/spray. **P273** Avoid release to the environment. **P280** Wear protective gloves/protective clothing/eye protection/face protection. **P305+P351+P338** IF IN EYES: Rinse cautiously with water for several minutes. Remove contact lenses, if present and easy to do. Continue rinsing. **P333+P313** If skin irritation or rash occurs: Get medical advice/attention. **P501** Dispose of contents/container in accordance with local/regional/national/international regulations.

Contains: 2-hydroxyethyl methacrylate, hydroxypropylmethacrylate, monoester with propane-1,2-diol, Diphenyl(2,4,6-trimethylbenzoyl) phosphin oxide, pigment paste

Lot number / Expiry date

The lot number and the expiry date are indicated on each SilaPrint model2 packaging. In case of claims please always indicate the lot number of the product. Do not use the product after the expiry date.


Inhalt/Contents: 1,0 kg

REF

267300

Achtung
Warning

Date of information: 02/2019

Printing date: 02/2019

SILADENT Dr. Böhme & Schöps GmbH

Im Klei 26

DE-38644 Goslar

Tel.: +49 (0) 53 21 / 37 79 - 0

Fax: +49 (0) 53 21 / 38 96 32

info@siladent.de

www.siladent.de