SilaPrint tray 385 nm/405nm

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Produktbeschreibung

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
Instructions for use

Product description

SilaPrint tray is a light-curing resin for the production of dental impression trays using 385 nm / 405 nm LED based stereolithography. For obligatory post curing of SilaPrint tray built parts, the high performance light curing unit Otoflash G171 is recommended, which is optimal balanced to the production process.

Processing

- During the application of SilaPrint tray we recommend wearing personal protective equipment f. i. suitable gloves, safety glasses etc.
- SilaPrint tray should be well shaken approx. 1 hour before use.
- Pour SilaPrint tray carefully into the pre-determined container of the production unit.
- Remove all bubbles with a cleaned object.
- Select the product specific parameter settings for your process.
- When building process is finished a direct post treatment is recommended.
- After the platform is taken up a dripping off time of approx. 10 minutes is recommended.
- Resin coated parts get clean with isopropanol (97 %) within a time of about 10 minutes. For a particular effective cleaning with isopropanol the use of explosion-protected ultrasonic units is recommended.
- Post-curing
Otoflash G171 with 2 x 3000 flashes under inert gas conditions (nitrogen), rotate components in between
- Impurity due to the construction or a break of the material and impurity due to operation mistakes cannot be excluded. Thanks to the low viscosity however it is possible to filtrate SilaPrint tray. It is recommended to regularly take out the container of the production unit, to homogenize and filtrate the content. Upcoming entrained bubbles get eliminated through a holding time of approx. 1 hour with SilaPrint tray.

To avoid damage of material quality does not expose the liquid material to light under any circumstances. Variations from the described manufacturing process may lead to modified mechanical characteristics and colour variations of the SilaPrint tray material.

SILADENT Dr. Böhme & Schöps GmbH

Im Klei 26
DE-38644 Goslar

Tel.: +49 (0) 53 21 / 37 79 - 0
Fax: +49 (0) 53 21 / 38 96 32

info@siladent.de
www.siladent.de

SilaPrint tray 385 nm/405 nm

Instructions for use

Possible dangers

Please note the safety instructions for processing SilaPrint tray:

Hazard statements: **H318** Causes serious eye damage. **H317** May cause an allergic skin reaction. **H412** Harmful to aquatic life with long lasting effects.

Precautionary statements: **P261** Avoid breathing dust/fume/gas/mist/vapours/spray. **P273** Avoid release to the environment. **P280** Wear protective gloves/eye protection/face protection. **P305+P351+P338** IF IN EYES: Rinse cautiously with water for several minutes. Remove contact lenses, if present and easy to do. Continue rinsing. **P310** Immediately call a POISON CENTER/ doctor. **P501** Dispose of contents/container in accordance with local/regional/national/international regulations.

Contains: 2-Propenoic acid, reaction product with Pentaerythrite, 1,4 Butanediol dimethacrylate, phenyl bis(2,4,6-trimethylbenzoyl)-phosphineoxide, Pigments

Lot number / Expiry date


The lot number and the expiry date are indicated on each SilaPrint tray packaging. In case of claims please always indicate the lot number of the product. Do not use the product after the expiry date.

Inhalt/Contents: 1,0 kg

 267304

Achtung
Warning

Date of information: 02/2019

Printing date: 02/2019

Hersteller / Manufacturer:

Dreve Dentamid GmbH
Max-Planck-Straße 31
59423 Unna / Germany

SILADENT Dr. Böhme & Schöps GmbH

Im Klei 26

DE-38644 Goslar

Tel.: +49 (0) 53 21 / 37 79 - 0

Fax: +49 (0) 53 21 / 38 96 32

info@siladent.de

www.siladent.de