# Splint Plus BioStar

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## Splint Plus BioStar

### Instructions for use

**Material composition:** Highly meshed (100 %) polycarbonate (PC)

**Indications:** Splint Plus BioStar dental milling blanks out of polycarbonate for the manufacture of splints, therapeutic splints, bite governors and drilling templates for short-term application in the oral cavity up to 30 days.

#### Contraindications

- intolerance to the contained components
- bridges with two molar pontics
- too little connector cross section
- model casting with clasps, bars, narrow sublingual bars and transversal connectors

<b>Physical characteristics:</b>	Flexural strength	(MPa)	≥ 98
	Flexural modulus	(MPa)	2.300
	Water absorption	( $\mu\text{g}/\text{mm}^3$ )	3.6
	Solubility	( $\mu\text{g}/\text{mm}^3$ )	0.2
	Residual monomer content	(%)	≤ 0.0015

**Safety information:** Please pay attention to the information in the instruction manual in its current version. Avoid the inhalation of dust particles during processing. Wear protection gloves, safety glasses and a dust mask to prevent irritations of eyes, skin and respiratory system.

**Handling and storage:** Only use clean, defect-free material out of undamaged packaging and make sure that the blanks are always stored in their original container in a dry and cool place (5-50°C), protected from light. Avoid contact with fluids.

**Processing:** Handling this medical device should only be allowed to trained staff. In general all constructions should be designed in a way to enable the compensation and the transfer of the mastication forces to the abutment teeth.

**Milling and further processing:**  Only use one-edged milling tools in the milling machine or wet grinding procedures!

**Further processing:** Tungsten carbide burs like the one you use for soft reline materials are suitable for rough adaption. For the finalization common cross-toothed millers can be used. Polishing should be carried out with regard to a low contact pressure, otherwise the resulting heat might lead to smearing or melting.

#### Cleaning

- ultrasonic bath (max. 40 °C) with mild cleaning agent
- process of sterilisation in an autoclave up to 121 °C

#### Do not use the following methods for cleaning:

- steam blasting
- organic solvents and acids (e.g. acetone, hydrochloric acid)
- highly concentrated alcohols
- chemical denture cleaners
- abrasive cleaning agents

**Mounting:** Bonding extra-orally: Only use special adhesives, adapted to the material you want to combine with the polycarbonate. Pay attention to the instructions of the adhesive manufacturer.

**Cementation:** For temporary cementations temporary cements and for permanent cementations self-adhesive glass ionomer cements are recommended.

**Disposal:** Remains of the Splint Plus BioStar blank can be disposed of in the regular household garbage.

**Final remarks:** All details are based on our own experiences and of details out of the literature which we have listed in the clinical assessment of the SILADENT. The statements and recommendations made do not lay any claim to completeness. We reserve the right to change these instructions of use.

**Further questions:** Please contact SILADENT technical services (Tel.: +49 (0) 53 21 – 37 79 25 / 26) or our sales representatives.

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