Splint Plus BioStar

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Indications:

Splint Plus BioStar are transparent dental milling blanks out of polycarbonate for the manufacture of splints, therapeutic splints, bite regulators and drilling templates for short-term application in the oral cavity up to 30 days.

Contraindications:

- intolerance to the contained components
- bridges with two molar pontics
- too little connector cross sections
- model casting with clasps, bars, narrow sublingual bars and transversal connectors

Safety information:

Please pay attention to the information in the instruction manual in its current version. Avoid the inhalation of dust particles during processing. Wear protection gloves, safety glasses and a dust mask to prevent irritations of eyes, skin and respiratory system.

Handling and storage:

Only use clean, defect-free material out of undamaged packaging and make sure that the blanks are always stored in their original container in a dry and cool place, protected from light. Avoid contact with fluids.

Milling and further processing:

Only use one-edged milling tools in the milling machine or wet grinding procedures!

Further processing:

Tungsten carbide burs like the one you use for soft reline materials are suitable for rough adaption. For the finalization common cross-toothed millers can be used. Polishing should be carried out with regard to a low contact pressure, otherwise the resulting heat might lead to smearing or melting.

Cleaning:

- ultrasonic bath (max. 40 °C) with mild cleaning agent
- process of sterilisation in an autoclave up to 121 °C

Do not use the following methods for cleaning:

- steam blasting
- organic solvents and acids (e.g. acetone, hydrochloric acid)
- highly concentrated alcohols
- chemical denture cleaners
- abrasive cleaning agents

Disposal:

Remains of the polycarbonate blank can be disposed of in the regular household garbage.

We reserve the right to make changes because of the continuous development of our products. Please find the current version of the instruction manual on our website: www.siladent.de This version replaces all previous versions.

Date of information: 11/2020

Printing date: 19.11.2020

Explanation of symbols

	Manufacturer		Date of manufacture
	Consult instructions for use		Use by
	Batch code		Catalogue No.
	Store dry		Protect from direct sun light
	Caution: U.S. Federal law restricts this device to sale by or on the order of a dentist.		
	Content (Quantity)		

Further questions:

Contact SILADENT technical services
(Tel.: +49 (0) 53 21 – 37 79 26)
or our sales representatives.

REV: 11/2020