## V-Alloy FG

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## V-Alloy FG

### Instructions for Use

**Application: Non precious alloy based on cobalt for removable dentures acc. DIN EN ISO 22674, type 5**

**Composition w<sub>i</sub>:**

Co 63.0% · Cr 30.0% · Mo 5.0%

**Other constituents:** Si, Mn, C

**Technic data:**

Proof stress 0.2 % (MPa)	745
Elongation at rupture (%)	13
Modulus of elasticity (GPa)	200
Vickers hardness VH 10	390
Density (g/cm <sup>3</sup> )	8.3
Melting range (°C)	1.300-1.370
Casting temperature (°C)	1.510

Tolerances in the composition of the alloy (% by mass) are within the range of current DIN standards.

**Investing and casting:**

V-Alloy FG can be used with phosphatebonded CrCo investments such as Granisit®, Micro, Jet 2000 or Granisit® XF Speed. Prepare the moulds for casting according to the SILADENT-System. Preheating range is 800 °C - 1000 °C. SILADENT recommends 920 °C. Heat soak at casting temperature for 40 minutes. Place the hot muffle and crucible only on a heat resistant underlayer (Risk of fire! Take care for safety distance to fire-endangered materials). Inductive cast: Glaze the crucible once before first use and heat it up in the preheating furnace. Afterwards a new glazing is not necessary. Start the cast when all the casting cubes have melted together and a final shadow runs over the molten metal just before the oxide layer disintegrates. When melting with a vacuum, the oxide layer in the vacuum disintegrates more quickly. Considerable experience is required when melting alloys automatically, and it is most important that the heating is set exactly according to the manufacturer's instructions. Accurate heat setting is essential to protect the properties of the alloy. Inaccurate heat setting impairs the alloy. Always use the reducing zone of the flame to melt the alloy. Ensure that casting is carried out as quickly as possible and use only preheated and cleaned crucibles. Overheating or a wrong preheating time produces poor results and may have a negative influence to the chemical and technical data.

**Prepare and polishing:**

To reduce the dust, soak the muffle in water before divesting. Do not knock against the casted model or against the casting cone during divesting. Risk of distortion! After sandblasting with aluminium oxide 250 m (distance to the blast nozzle 10 – 15 mm), the framework can be prepared with rotating stones, milling tools, with rubber abrasives, brushes, buffing wheels and polishing paste.

**Soldering and welding:**

Grind the affected parts of the corrosion resistant alloy with cross-cut carbide burs (do not sandblast). Heat up the soldering object to the temperature of the solder (e.g. SILADENT CoCr special-solder). If the necessary temperature is reached, the solder, while using soldering flux, should flow into the soldered joint. If the solder flows in small bowls the temperature is not high enough. For lasering we recommend SILADENT CrCo laser welding rods.

**Safety hint:**

Metal dust is hazardous to health. For finishing and sandblasting use a suitable extraction system and / or face mask.

**Contraindications, side effects and interactions with other dental alloys:**

Incompatibility reactions are very seldom, when working instructions are strictly considered. In cases of hypersensitivity (allergy) to the constituents of the alloy, discontinue its use. In individual cases, electrochemically induced local dysaesthesia has been reported. Galvanic effects can occur under proximal or antagonistic contact with dentures of different alloys.

**Storage conditions:**

No special storage conditions are required.

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**Manufacturer:** The Argen Corporation · 5855 Oberlin Drive, San Diego · California 92121, USA

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Date of information: 10/2020

Printing date: 10/2020

**Further questions:** Contact SILADENT technical services  
(Tel.: +49 (0) 53 21 – 37 79 19)  
or our sales representatives.

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**Distribuidor:** SILADENT Dr. Böhme & Schöps GmbH, Im Klei 26, DE 38644 Goslar
**Fabricante:** The Argen Corporation · 5855 Oberlin Drive, San Diego · California 92121, USA
The Argen Corporation rotula el envase comercial con **CE 2797**
**SILADENT Dr. Böhme & Schöps GmbH**