

## Einweg-Pressstempel

DIN EN ISO 15912

**Liefergebände:** REF 102650 Packung á 50 Stück, Ø 12 mm  
REF 102651 Packung á 50 Stück, Ø 13 mm

### Vorteile:

- kann ohne Vorwärmen in die vorgeheizte Muffel gegeben werden
- einfaches Handling
- erspart ein zeitraubendes Abstrahlen von wiederverwendbaren Pressstempeln
- keine Mikrosprünge
- geeignet für alle Presskeramiken
- in 2 Größen lieferbar: Ø 12 mm (z. B. Empress®/Ivoclar) und Ø 13 mm (z.B. e.max® / Ivoclar)

### Verarbeitungshinweise:

Gemäß den Verarbeitungsempfehlungen der Presskeramiklieferanten verarbeiten. Nicht Vorwärmen. Markierungspunkt oben!

Achtung: Vor Feuchtigkeit schützen!

Stand der Information: 11/2014

## Disposable Plunger

DIN EN ISO 15912

**Packing sizes:** REF 102650 carton á 50 pieces, Ø 12 mm  
REF 102651 carton á 50 pieces, Ø 13 mm

### Advantages:

- No preheating of the plunger!
- Easy to handle
- no more time-consuming divesting of reusable plungers
- no micro cracks in the sprues
- suitable for all pressable ceramic systems
- available in 2 diameters: Ø 12 mm (e.g. Empress®/Ivoclar) and Ø 13 mm (e. g. e.max®/Ivoclar).

### Instructions for use:

Follow the instruction of the manufacturer for the used pressable ceramics. Do not preheat the plungers! Marker on the top.

Caution: Protect from humidity!

Date of information: 11/2014