

According to 1907/2006/EG, Article 31

Revision: 31.08.2017

Version number 3

Page 1 of 7

Printing date: 06.12.2022

SilaPrint cast**SECTION 1: Identification of the substance/mixture and of the company/undertaking** ***1.1 Product identifier****Commercial product name:** SilaPrint cast**1.2 Relevant identified uses of the substance or mixture and uses advised against:** No further relevant information available.**Application of the substance / the mixture:** Material based on methacrylate resin for DLP systems with 385nm / 405nm LED for the manufacturing of profiles that burn out without any residue for the dental casting technology**1.3 Details of the supplier of the safety data sheet**

Manufacturer/Supplier: SILADENT Dr. Böhme & Schöps GmbH
Street / mailbox: Im Klei 26
Country code. / postal code / city: DE - 38644 Goslar
Phone: +49 (0) 53 21 / 37 79 - 0
Fax: +49 (0) 53 21 / 38 96 32
E-mail / Website: info@siladent.de / www.siladent.de
Further information obtainable from: SILADENT Dr. Böhme & Schöps GmbH

1.4 Emergency telephone number**SILADENT Dr. Böhme & Schöps GmbH:** +49 (0) 53 21 / 37 79 - 0 (Mon-Fri. 8 a.m. - 4 p.m.)**SECTION 2: Hazards identification****2.1 Classification of the substance or mixture****Classification according to Regulation (EC) No 1272/2008:**

GHS07

Skin Irrit. 2	H315	Causes skin irritation.
Eye Irrit. 2	H319	Causes serious eye irritation.
Skin Sens. 1	H317	May cause an allergic skin reaction.
STOT SE 3	H335	May cause respiratory irritation.

2.2 Label elements:**Labelling according to Regulation (EC) No 1272/2008:**

The product is classified and labelled according to the CLP regulation.

Hazard pictograms:

GHS07

Signal word:

Warning

Hazard-determining components of labelling:Diurethandimethacrylate
triethylene glycol dimethacrylate
phenyl bis(2,4,6-trimethylbenzoyl)-phosphineoxide**Hazard statements:**

H315	Causes skin irritation.
H319	Causes serious eye irritation.
H317	May cause an allergic skin reaction.
H335	May cause respiratory irritation.

Precautionary statements:

According to 1907/2006/EG, Article 31

Revision: 31.08.2017

Version number 3

Page 2 of 7

Printing date: 06.12.2022

SilaPrint cast

P261	Avoid breathing dust/fume/gas/mist/vapours/spray.
P280:	Wear protective gloves / eye protection.
P304+P340	IF INHALED: Remove person to fresh air and keep comfortable for breathing.
P305+P351+P338:	IF IN EYES: Rinse cautiously with water for several minutes. Remove contact lenses, if present and easy to do. Continue rinsing.
P405	Store locked up.
P501	Dispose of contents/container in accordance with local/regional/national/international regulations

2.3 Other hazards:**Results of PBT and vPvB assessment:**

PBT:	Not applicable.
vPvB:	Not applicable.

SECTION 3: Composition/information on ingredients**3.1 Chemical characterisation: Mixtures**

Description: Mixture of substances listed below with nonhazardous additions.

Dangerous components:		
CAS: 72869-86-4 EINECS: 276-957-5	Diurethandimethacrylate	50-100 %
	⚠ Skin Sens. 1, H317	
CAS: 109-16-0	triethylene glycol dimethacrylate	25-50%
	⚠ Skin Irrit. 2, H315; Eye Irrit. 2, H319; Skin Sens. 1, H317; STOT SE 3, H335	
CAS: 162881-26-7 ELINCS: 423-340-5	phenyl bis(2,4,6-trimethylbenzoyl)-phosphineoxide	≤ 2.5%
	⚠ Skin Sens. 1, H317; Aquatic Chronic 4, H413	

Additional information:

For the wording of the listed risk phrases refer to section 16.

SECTION 4: First aid measures**4.1 Description of first aid measures**

General information:	Immediately remove any clothing soiled by the product.
After inhalation:	Supply fresh air and to be sure call for a doctor. In case of unconsciousness place patient stably in side position for transportation.
After skin contact:	Immediately wash with water and soap and rinse thoroughly. If skin irritation continues, consult a doctor.
After eye contact:	Rinse opened eye for several minutes under running water. Then consult a doctor.
After swallowing:	Do not induce vomiting. Call a doctor immediately.

4.2 Most important symptoms and effects, both acute and delayed:

No further relevant information available.

4.3 Indication of any immediate medical attention and special treatment needed:

No further relevant information available.

SECTION 5: Firefighting measures**5.1 Extinguishing media**

According to 1907/2006/EG, Article 31

Revision: 31.08.2017

Version number 3

Page 3 of 7

Printing date: 06.12.2022

SilaPrint cast

Suitable extinguishing agents: CO₂, powder or water spray. Fight larger fires with water spray or alcohol resistant foam.

For safety reasons unsuitable extinguishing agents: Water with full jet

5.2 Special hazards arising from the substance or mixture: No further relevant information available.

**5.3 Advice for firefighters:
Protective equipment:** No special measures required.

SECTION 6: Accidental release measures

6.1 Personal precautions, protective equipment and emergency procedures: Avoid contact with skin and eyes.

6.2 Environmental precautions: Do not allow product to reach sewage system or any water course.

6.3 Methods and material for containment and cleaning up: Absorb with liquid-binding material (sand, diatomite, acid binders, universal binders, sawdust).
Use neutralising agent.
Dispose contaminated material as waste according to item 13.
Ensure adequate ventilation.

6.4 Reference to other sections See Section 7 for information on safe handling.
See Section 8 for information on personal protection equipment. See Section 13 for disposal information.

SECTION 7: Handling and storage

7.1 Precautions for safe handling: Avoid contact with skin and eyes.
Ensure good ventilation/exhaustion at the workplace.
Prevent formation of aerosols.

Information about fire - and explosion protection: No special measures required.

**7.2 Conditions for safe storage, including any incompatibilities
Storage:**

Requirements to be met by storerooms and receptacles: Protect from heat and direct sunlight.

Information about storage in one common storage facility: Not required.

Further information about storage conditions: Keep container tightly sealed.

7.3 Specific end use(s): No further relevant information available.

SECTION 8: Exposure controls/personal protection

Additional information about design of technical facilities: No further data; see item 7.

8.1 Control parameters:

SilaPrint cast**Ingredients with limit values that require monitoring at the workplace:**

The product does not contain any relevant quantities of materials with critical values that have to be monitored at the workplace.

8.2 Additional information:**Exposure controls****Personal protective equipment:****General protective and hygienic measures:**

The lists valid during the making were used as basis.

Keep away from foodstuffs, beverages and feed.
Immediately remove all soiled and contaminated clothing
Wash hands before breaks and at the end of work.
Avoid contact with the eyes and skin.

Respiratory protection:

In case of brief exposure or low pollution use respiratory filter device. In case of intensive or longer exposure use self-contained respiratory protective device.

Protection of hands:

Protective gloves

Protective gloves should be changed regularly, especially after intensive contact with the product.

For every workplace a suitable type of protective gloves must be selected. The glove material has to be impermeable and resistant to the product/ the substance/ the preparation.

Material of gloves:

Nitrile rubber, NBR (0,4 mm)

As there are many different conditions in every day work these indications can only serve as an aid to orientation for the selection of suitable gloves for the handling of chemical products. By no means can they replace qualifying examinations by the end-user. These recommendations only apply to the product mentioned in the safety data sheet. When mixing with other substances or under conditions deviant from norm EN 374 a manufacturer of CE-approved gloves should be referred to. The selection of the suitable gloves does not only depend on the material, but also on further marks of quality and varies from manufacturer to manufacturer. As the product is a preparation of several substances, the resistance of the glove material cannot be calculated in advance and has therefore to be checked prior to the application.

Penetration time of glove material:

The exact break through time has to be found out by the manufacturer of the protective gloves and has to be observed. The determined penetration times according to EN 374 part III are not performed under practical conditions. Therefore a maximum wearing time, which corresponds to 50% of the penetration time, is recommended.

Eye protection:

Tightly sealed goggles

Body protection:

Protective work clothing

According to 1907/2006/EG, Article 31

Revision: 31.08.2017

Version number 3

Page 5 of 7

Printing date: 06.12.2022

SilaPrint cast**SECTION 9: Physical and chemical properties****9.1 Information on basic physical and chemical properties****General Information****Appearance:****Form:**

Fluid

Colour:

Reddish transparent

Odour:

Characteristic

Odour threshold:

Not determined.

pH-value:

Not determined.

Change in condition**Melting point/Melting range:**

Undetermined.

Boiling point/Boiling range:

201 °C

Flash point:

151 °C

Flammability (solid, gas):

Not applicable.

Decomposition temperature:

445 °C

Auto-ignition temperature:

Product is not selfigniting.

Explosive properties:

Product does not present an explosion hazard.

Explosion limits:**Lower:**

Not determined.

Upper:

Not determined.

Oxidising properties

Not determined.

Vapour pressure:

Not determined.

Density at 20 °C:1.1 g/cm³**Relative density:**

Not determined.

Vapour density:

Not determined.

Evaporation rate:

Not determined.

Solubility in / Miscibility with water:

Not miscible or difficult to mix.

Partition coefficient: n-octanol/water:

Not determined.

Viscosity:**Dynamic at 23 °C:**

200 mPas

Kinematic:

Not determined.

Solvent content:**Organic solvents:**

0.0 %

9.2 Other information:

No further relevant information available.

SECTION 10: Stability and reactivity**10.1 Reactivity:**

No dangerous reactions if used according to specifications.

10.2 Chemical stability:

Stable if used according to specifications.

Thermal decomposition / conditions to be avoided:

Protect from heat and direct sunlight.

No decomposition if used according to specifications.

10.3 Possibility of hazardous reactions:

No dangerous reactions known.

10.4 Conditions to avoid:

No further relevant information available.

10.5 Incompatible materials:

No further relevant information available.

10.6 Hazardous decomposition products:

No dangerous decomposition products known.

SECTION 11: Toxicological information**11.1 Information on toxicological effects****Acute toxicity:**

Based on available data, the classification criteria are not met.

Primary irritant effect:

According to 1907/2006/EG, Article 31

Revision: 31.08.2017

Version number 3

Page 6 of 7

Printing date: 06.12.2022

SilaPrint cast

Skin corrosion/irritation:	Causes skin irritation.
Serious eye damage/irritation:	Causes serious eye irritation.
Respiratory or skin sensitisation:	May cause an allergic skin reaction.
CMR effects (carcinogenicity, mutagenicity and toxicity for reproduction)	
Germ cell mutagenicity:	Based on available data, the classification criteria are not met.
Carcinogenicity:	Based on available data, the classification criteria are not met.
Reproductive toxicity:	Based on available data, the classification criteria are not met.
STOT-single exposure:	May cause respiratory irritation.
STOT-repeated exposure:	Based on available data, the classification criteria are not met.
Aspiration hazard:	Based on available data, the classification criteria are not met.

SECTION 12: Ecological information

12.1 Toxicity	
Aquatic toxicity:	No further relevant information available.
12.2 Persistence and degradability:	No further relevant information available.
12.3 Bioaccumulative potential:	No further relevant information available.
12.4 Mobility in soil:	No further relevant information available.
Additional ecological information:	
General notes:	Do not allow product to reach ground water, water course or sewage system.
12.5 Results of PBT and vPvB assessment:	
PBT:	Not applicable.
vPvB:	Not applicable.
12.6 Other adverse effects:	No further relevant information available.

SECTION 13: Disposal considerations

13.1 Waste treatment methods	
Recommendation:	Must not be disposed together with household garbage. Do not allow product to reach sewage system. Disposal must be made according to official regulations.

European waste catalogue	
HP 4	Irritant - skin irritation and eye damage
HP 5	Specific Target Organ Toxicity (STOT)/Aspiration Toxicity
HP 13	Sensitising

13.2 Uncleaned packaging:	
Recommendation:	Disposal must be made according to official regulations.

SECTION 14: Transport information

14.1 UN-Number	
ADR, ADN, IMDG, LATA	Void
14.2 UN proper shipping name	
ADR, AND, IMDG, LATA	Void
14.3 Transport hazard class(es)	
ADR, ADN, IMDG, LATA	
Class	Void
14.4 Packing group	

According to 1907/2006/EG, Article 31

Revision: 31.08.2017

Version number 3

Page 7 of 7

Printing date: 06.12.2022

SilaPrint cast**ADR, IMDG, LATA**

Void

14.5 Environmental hazards:**Marine pollutant:**

No

14.6 Special precautions for user

Not applicable

14.7 Transport in bulk according to Annex II of MARPOL73/78 and the IBC Code UN "Model Regulation":

Not applicable.

Void

SECTION 15: Regulatory information**15.1 Safety, health and environmental regulations/legislation specific for the substance or mixture Directive 2012/18/EU****Named dangerous substances - ANNEX I:**

None of the ingredients is listed.

REGULATION (EC) No 1907/2006 ANNEX XVII:

Conditions of restriction: 3

15.2 Chemical safety assessment:

A Chemical Safety Assessment has not been carried out.

SECTION 16: Other information

This information is based on our present knowledge. However, this shall not constitute a guarantee for any specific product features and shall not establish a legally valid contractual relationship.

16.1 Relevant phrases:

H315 Causes skin irritation.
 H317 May cause an allergic skin reaction.
 H319 Causes serious eye irritation.
 H335 May cause respiratory irritation.
 H413 May cause long lasting harmful effects to aquatic life.

16.2 Abbreviations and acronyms:

RID: Règlement international concernant le transport des marchandises dangereuses par chemin de fer (Regulations Concerning the International Transport of Dangerous Goods by Rail)
 ICAO: International Civil Aviation Organisation
 ADR: Accord européen sur le transport des marchandises dangereuses par Route (European Agreement concerning the International Carriage of Dangerous Goods by Road)
 IMDG: International Maritime Code for Dangerous Goods
 IATA: International Air Transport Association
 GHS: Globally Harmonised System of Classification and Labelling of Chemicals
 EINECS: European Inventory of Existing Commercial Chemical Substances
 ELINCS: European List of Notified Chemical Substances
 CAS: Chemical Abstracts Service (division of the American Chemical Society)
 PBT: Persistent, Bioaccumulative and Toxic
 vPvB: very Persistent and very Bioaccumulative
 Skin Irrit. 2: Skin corrosion/irritation – Category 2
 Eye Irrit. 2: Serious eye damage/eye irritation – Category 2
 Skin Sens. 1: Skin sensitisation – Category 1
 Aquatic Chronic 4: Hazardous to the aquatic environment - long-term aquatic hazard – Category 4

* Data compared to the previous version altered.