

Splint Plus BioStar**SECTION 1: Identification of the substance/mixture and of the company/undertaking**

1.1 Product identifier:**Commercial product name:** Splint Plus BioStar**1.2 Relevant identified uses of the substance or mixture and uses advised against****Use of the substance/mixture:**

Splint Plus BioStar are transparent dental milling blanks out of polycarbonate for the manufacture of splints, therapeutic splints, bite regulators and drilling templates for short-term application in the oral cavity up to 30 days.

Uses advised against:

No information available.

1.3 Details of the supplier of the safety data sheet**Manufacturer/Supplier:**

SILADENT Dr. Böhme & Schöps GmbH

Street / mailbox:

Im Klei 26

Country code. / postal code / city:

DE - 38644 Goslar

Phone:

+49 (0) 53 21 / 37 79 - 0

Fax:

+49 (0) 53 21 / 38 96 32

E-mail / Website:info@siladent.de / www.siladent.de**1.4 Emergency telephone number**

SILADENT Dr. Böhme & Schöps GmbH: +49 (0) 53 21 / 37 79 - 0 (Mon-Fri. 8 a.m. - 4 p.m.)

SECTION 2: Hazards identification

2.1 Classification of the substance or mixture**Regulation (EC) No. 1272/2008:**

This substance is not classified as hazardous according to Regulation (EC) No. 1272/2008.

2.2 Label Elements**2.3 Other hazards:**

No information available.

SECTION 3: Composition/information on ingredients

3.1 Substances**Chemical characterization:**

No information available.

Further Information:

No information available

SECTION 4: First aid measures

4.1 Description of first aid measures**After inhalation:**

Remove person to fresh air and keep comfortable for breathing. In case of respiratory tract irritation, consult a physician.

After contact with skin:

After contact with skin, wash immediately with plenty of water and soap. If skin irritation occurs: Get medical advice/attention.

After contact with eyes:

Rinse immediately carefully and thoroughly with eye-bath or water. In case of eye irritation consult an ophthalmologist.

After ingestion:

IF SWALLOWED: Call a doctor if you feel unwell.

4.2 Most important symptoms and effects, both acute and delayed:

No information available.

4.3 Indication of any immediate medical attention and special treatment needed:

No information available.

Splint Plus BioStar

SECTION 5: Firefighting measures

- 5.1 Extinguishing media**
- Suitable extinguishing agent:** Dry extinguishing powder. Carbon dioxide (CO₂). Alcohol resistant foam. Water spray jet.
- Unsuitable extinguishing media:** High power water jet.
- 5.2 Special hazards arising from the substance or mixture:** In case of fire may be liberated: Carbon monoxide, Carbon dioxide, Nitrogen oxides (NO_x). Possible in traces: Hydrogen cyanide (hydrocyanic acid) Do not inhale explosion and combustion gases.
- 5.3 Advice for firefighters:** Co-ordinate fire-fighting measures to the fire surroundings. In case of fire: Wear self-contained breathing apparatus.
- Additional information:** Do not allow run-off from fire-fighting to enter drains or water courses.

SECTION 6: Accidental release measures

- 6.1 Personal precautions, protective equipment and emergency procedures:** Use personal protection equipment.
- 6.2 Environmental precautions:** Do not allow to enter into surface water or drains.
- 6.3 Methods and material for containment and cleaning up:** Take up mechanically. Avoid dust formation. Collect in closed and suitable containers for disposal.
- 6.4 Reference to other sections:** See protective measures under point 7 and 8. Disposal: see section 13.

SECTION 7: Handling and storage

- 7.1 Precautions for safe handling**
- Advice on safe handling:** See section 8.
- Advice on protection against fire and explosion:** No special fire protection measures are necessary.
- Further information on handling:** Avoid dust formation. Protect from moisture. Provide adequate ventilation as well as local exhaust at critical locations. To follow: Occupational exposure limit values.
- 7.2 Conditions for safe storage, including any incompatibilities**
- Requirements for storage rooms and vessels:** Keep/Store only in original container.
- Advice on storage compatibility:** Keep the packing dry and well-sealed to prevent contamination and absorption of humidity.
- Further information on storage conditions:** Keep away from: Frost, Heat Handle with care - avoid bumps, friction and impact.
- 7.3 Specific end use(s):** Splint Plus BioStar are transparent dental milling blanks out of polycarbonate for the manufacture of splints, therapeutic splints, bite regulators and drilling templates for short-term application in the oral cavity up to 30 days.

Splint Plus BioStar

SECTION 8: Exposure controls/personal protection

8.1 Control parameters

8.2 Exposure controls

Appropriate engineering controls:	Technical ventilation of workplace.
Protective and hygiene measures:	Work in well-ventilated zones or use proper respiratory protection. Only wear fitting, comfortable and clean protective clothing. Wash hands before breaks and after work. Separate storage of work clothes.
Eye/face protection:	Suitable eye protection: Eye glasses with side protection, goggles
Hand protection:	Suitable gloves type: NBR (Nitrile rubber) DIN EN 374, Butyl caoutchouc (butyl rubber) DIN EN 374, PVC (Polyvinyl chloride) Thickness of the glove material: >0,5 mm Wear cotton under mitten if possible. Breakthrough times and swelling properties of the material must be taken into consideration.
Skin protection:	Protective clothing.
Respiratory protection:	Particle filter device (DIN EN 143)

SECTION 9: Physical and chemical properties

9.1 Information on basic physical and chemical properties

Physical state:	solid	
Colour:	transparent	
Odour:	odourless	
pH-Value:	not determined	Test method
Changes in the physical state		
Melting point:	not determined	
Initial boiling point and boiling range:	not determined	
Sublimation point:	not determined	
Softening point:	not determined	
Pour point:	not determined	
Flash point:	not determined	
Flammability		
Solid:	>450 °C	
Gas:	not applicable	
Explosive properties:		
Lower explosion limits:	No information available.	
Upper explosion limits:	not applicable	
Auto-ignition temperature		
Solid:	not determined	
Gas:	not determined	
Decomposition temperature:	not determined	
Oxidizing properties:		
Vapour pressure:	not determined	
Density:	not determined	
Bulk density:	600-700 kg/m ³	
Water solubility:	insoluble	
Solubility in other solvents:	No information available.	

Splint Plus BioStar

Partition coefficient:	not applicable
Viscosity/dynamic:	not determined
Vapour density:	not determined
Evaporation rate:	not determined

9.2 Other information: No information available.

SECTION 10: Stability and reactivity

10.1 Reactivity: No information available.

10.2 Chemical stability: No information available.

10.3 Possibility of hazardous reactions: In case of fire may be liberated: Gases/vapours, harmful.

10.4 Conditions to avoid: No information available.

10.5 Incompatible materials: No information available.

10.6 Hazardous decomposition products: Carbon dioxide (CO₂), CO

SECTION 11: Toxicological information

11.1 Information on toxicological effects
Toxicocinetics, metabolism and distribution: No information available.

Acute toxicity: Based on available data, the classification criteria are not met.

Irritation and corrosivity: Based on available data, the classification criteria are not met.

Sensitising effects: Based on available data, the classification criteria are not met.

Carcinogenic/mutagenic/toxic effects for reproduction: Based on available data, the classification criteria are not met.

STOT-single exposure: Based on available data, the classification criteria are not met.

STOT-repeated exposure: Based on available data, the classification criteria are not met.

Aspiration hazard: Based on available data, the classification criteria are not met.

SECTION 12: Ecological information

12.1 Toxicity: No information available.

12.2 Persistence and degradability: The product is: Not readily biodegradable (according to OECD criteria).

12.3 Bioaccumulative potential: No information available.

12.4 Mobility in soil: No information available.

12.5 Results of PBT and vPvB assessment: The substances in the mixture do not meet the PBT/vPvB criteria according to REACH, annex XIII.

12.6 Other adverse effects: No information available.

Splint Plus BioStar**SECTION 13: Disposal considerations****13.1 Waste treatment methods****Advice on disposal:**

Dispose of waste according to applicable legislation.

Contaminated packaging:

Dispose of waste according to applicable legislation. Non-contaminated packages may be recycled.

SECTION 14: Transport information**Land transport (ADR/RID)****14.1 UN number:**

No dangerous good in sense of this transport regulation.

14.2 UN proper shipping name:

No dangerous good in sense of this transport regulation.

14.3 Transport hazard class(es):

No dangerous good in sense of this transport regulation.

14.4 Packing group:

No dangerous good in sense of this transport regulation.

Inland waterways transport (ADN)**14.1 UN number:**

No dangerous good in sense of this transport regulation.

14.2 UN proper shipping name:

No dangerous good in sense of this transport regulation.

14.3 Transport hazard class(es):

No dangerous good in sense of this transport regulation.

14.4 Packing group:

No dangerous good in sense of this transport regulation.

Marine transport (IMDG)**14.1 UN number:**

No dangerous good in sense of this transport regulation.

14.2 UN proper shipping name:

No dangerous good in sense of this transport regulation.

14.3 Transport hazard class(es):

No dangerous good in sense of this transport regulation.

14.4 Packing group:

No dangerous good in sense of this transport regulation.

Air transport (ICAO)**14.1 UN number:**

No dangerous good in sense of this transport regulation.

14.2 UN proper shipping name:

No dangerous good in sense of this transport regulation.

14.3 Transport hazard class(es):

No dangerous good in sense of this transport regulation.

14.4 Packing group:

No dangerous good in sense of this transport regulation.

14.5 Environmental hazards**ENVIRONMENTALLY HAZARDOUS:**

No

14.6 Special precautions for user:

No dangerous good in sense of this transport regulation.

14.7 Transport in bulk according to Annex II of Marpol and the IBC Code:

No dangerous good in sense of this transport regulation.

SECTION 15: Regulatory information**15.1 Safety, health and environmental regulations/legislation specific for the substance or mixture****National regulatory information****Water contaminating class (D):**

- - not water contaminating

15.2 Chemical safety assessment:

Chemical safety assessments for substances in this mixture were not carried out.

SECTION 16: Other information**Abbreviations and acronyms**

ADR: Accord européen sur le transport des marchandises dangereuses par Route (European Agreement concerning the International Carriage of Dangerous Goods by Road)

RID: Règlement international conernat le transport des marchandises dangereuses par chemin de fer (Regulations Concerning the International Transport of Dangerous Goods by Rail)

IMDG: International Maritime Code for Dangerous Goods

IATA: International Air Transport Association

IATA-DGR: Dangerous Goods Regulations by the "International Air Transport Association" (IATA)

Splint Plus BioStar

	ICAO: International Civil Aviation Organization
ICAO-TI:	Technical Instructions by the "International Civil Aviation Organization" (ICAO)
CAS:	Chemical Abstracts Service (division of the American Chemical Society)
GHS:	Globally Harmonized System of Classification and Labelling of Chemicals
CLP:	Regulation on Classification, Labelling and Packaging of Substances and Mixtures,
LC50:	Lethal concentration, 50 percent
LD50:	Lethal dose, 50 percent
EC50:	Effective concentration, 50 percent
DNEL:	Derived No Effect Level
PNEC:	Predicted No Effect Concentration
PBT:	Persistent, Bioaccumulative and Toxic
vPvB:	very Persistent and very Bioaccumulative

Further Information

The above information describes exclusively the safety requirements of the product and is based on our present-day knowledge. The information is intended to give you advice about the safe handling of the product named in this safety data sheet, for storage, processing, transport and disposal. The information cannot be transferred to other products. In the case of mixing the product with other products or in the case of processing, the information on this safety data sheet is not necessarily valid for the new made-up material