

Zirkon BioStar, Zirkon BioStar Colour**SECTION 1: Identification of the substance/mixture and of the company/undertaking**

1.1 Product identifier:**Commercial product name:** Zirkon BioStar, Zirkon BioStar Colour**1.2 Relevant identified uses of the substance or mixture and uses advised against****Use of the substance/mixture:** Zirkon BioStar, Zirkon BioStar Colour are dental milling blanks made of yttria-stabilized zirconia for the manufacture of dental restorations. They can be used for all dental constructions with up to two pontics in anterior and posterior region.**Uses advised against:** No information available.**1.3 Details of the supplier of the safety data sheet****Manufacturer:**

Company name: Dental Direkt GmbH
Street: Industriezentrum 106-108
Place: DE-32139 Spenge
Telephone: +49-(0)5225 - 8 63 19-0
Telefax: +49-(0)5225 - 8 63 19-99
e-mail:
Contact person: main office - Telephone: +49-(0)5225 - 8 63 19-0
Internet: www.dentaldirekt.de
Responsible Department: info@dentaldirekt.de

Supplier:

Street / mailbox: SILADENT Dr. Böhme & Schöps GmbH
Im Klei 26
Country code. / postal code / city: DE - 38644 Goslar
Phone: +49 (0) 53 21 / 37 79 - 0
Fax: +49 (0) 53 21 / 38 96 32
E-mail / Website: info@siladent.de / www.siladent.de

1.4 Emergency telephone number**Further Information:** +49 (0) 761 19240 (VIZ Freiburg)

No information available.

SECTION 2: Hazards identification

2.1 Classification of the substance or mixture**Regulation (EC) No. 1272/2008:** This mixture is not classified as hazardous according to Regulation (EC) No. 1272/2008.**2.2 Label Elements****2.3 Other hazards:** No information available.**SECTION 3: Composition/information on ingredients**

3.2 Mixtures**Chemical characterization:** No information available.**Further Information:** No information available**SECTION 4: First aid measures**

4.1 Description of first aid measures**After inhalation:** Remove person to fresh air and keep comfortable for breathing. In case of respiratory tract irritation, consult a physician.

Zirkon BioStar, Zirkon BioStar Colour

After contact with skin:	After contact with skin, wash immediately with plenty of water and soap. If skin irritation occurs: Get medical advice/attention.
After contact with eyes:	Rinse immediately carefully and thoroughly with eye-bath or water. In case of eye irritation consult an ophthalmologist.
After ingestion:	Rinse mouth immediately and drink plenty of water.
4.2 Most important symptoms and effects, both acute and delayed:	No information available.
4.3 Indication of any immediate medical attention and special treatment needed:	No information available.

SECTION 5: Firefighting measures

5.1 Extinguishing media	
Suitable extinguishing agent:	Dry extinguishing powder. Carbon dioxide (CO2). Alcohol resistant foam. Water spray jet.
Unsuitable extinguishing media:	High power water jet.
5.2 Special hazards arising from the substance or mixture:	No information available.
5.3 Advice for firefighters:	Co-ordinate fire-fighting measures to the fire surroundings. The product itself does not burn. In case of fire: Wear self-contained breathing apparatus.
Additional information:	Do not allow run-off from fire-fighting to enter drains or water courses.

SECTION 6: Accidental release measures

6.1 Personal precautions, protective equipment and emergency procedures:	Use personal protection equipment.
6.2 Environmental precautions:	Do not allow to enter into surface water or drains.
6.3 Methods and material for containment and cleaning up:	Take up mechanically. Avoid dust formation. Collect in closed and suitable containers for disposal.
6.4 Reference to other sections:	See protective measures under point 7 and 8. Disposal: see section 13.

SECTION 7: Handling and storage

7.1 Precautions for safe handling	
Advice on safe handling:	See section 8.
Advice on protection against fire and explosion:	No special fire protection measures are necessary.
Further information on handling:	Avoid dust formation. Protect from moisture. Provide adequate ventilation as well as local exhaustion at critical locations. To follow: Occupational exposure limit values.
7.2 Conditions for safe storage, including any incompatibilities	

Zirkon BioStar, Zirkon BioStar Colour

Requirements for storage rooms and vessels:	Keep/Store only in original container.
Advice on storage compatibility:	Keep the packing dry and well-sealed to prevent contamination and absorption of humidity.
Further information on storage conditions:	Keep away from: Frost, Heat Handle with care - avoid bumps, friction and impact.
7.3 Specific end use(s):	Observe instructions for use.

SECTION 8: Exposure controls/personal protection

8.1 Control parameters

8.2 Exposure controls

Appropriate engineering controls:	Technical ventilation of workplace.
Protective and hygiene measures:	Work in well-ventilated zones or use proper respiratory protection. Only wear fitting, comfortable and clean protective clothing. Wash hands before breaks and after work. Separate storage of work clothes.
Eye/face protection:	Suitable eye protection: Eye glasses with side protection, goggles
Hand protection:	Suitable gloves type: NBR (Nitrile rubber) DIN EN 374, Butyl caoutchouc (butyl rubber) DIN EN 374, NR (natural rubber, natural latex). Wear cotton under mitten if possible. Breakthrough times and swelling properties of the material must be taken into consideration.
Skin protection:	Protective clothing.
Respiratory protection:	Particle filter device (DIN EN 143)

SECTION 9: Physical and chemical properties

9.1 Information on basic physical and chemical properties

Physical state:	solid	
Colour:	white	
Odour:	odourless	
pH-Value:	not applicable	Test method
Changes in the physical state		
Melting point:	not determined	
Initial boiling point and boiling range:	not determined	
Sublimation point:	not determined	
Softening point:	not determined	
Pour point:	not determined	
Flash point:	not applicable	
Flammability		
Solid:	not applicable	
Gas:	not applicable	
Explosive properties:	not explosive according to EU A.14	
Lower explosion limits:	not applicable	

Zirkon BioStar, Zirkon BioStar Colour

Upper explosion limits:	not applicable
Auto-ignition temperature	
Solid:	not determined
Gas:	not determined
Decomposition temperature:	not determined
Oxidizing properties:	Not oxidizing.
Vapour pressure:	not determined
Density:	not determined
Bulk density:	1.000-1.500 kg/m ³
Water solubility:	insoluble
Solubility in other solvents:	No information available.
Partition coefficient:	not applicable
Viscosity/dynamic:	not determined
Vapour density:	not determined
Evaporation rate:	not determined

9.2 Other information: No information available.

SECTION 10: Stability and reactivity

10.1 Reactivity:	No information available.
10.2 Chemical stability:	No information available.
10.3 Possibility of hazardous reactions:	No information available.
10.4 Conditions to avoid:	No information available.
10.5 Incompatible materials:	No information available.
10.6 Hazardous decomposition products:	Not applicable.

SECTION 11: Toxicological information

11.1 Information on toxicological effects	
Toxicokinetics, metabolism and distribution:	No information available.
Acute toxicity:	Based on available data, the classification criteria are not met.
Irritation and corrosivity:	Based on available data, the classification criteria are not met.
Sensitising effects:	Based on available data, the classification criteria are not met.
Carcinogenic/mutagenic/toxic effects for reproduction:	Based on available data, the classification criteria are not met.
STOT-single exposure:	Based on available data, the classification criteria are not met.
STOT-repeated exposure:	Based on available data, the classification criteria are not met.
Aspiration hazard:	Based on available data, the classification criteria are not met.

SECTION 12: Ecological information

12.1 Toxicity:	No information available.
12.2 Persistence and degradability:	No information available.
12.3 Bioaccumulative potential:	No information available.

Zirkon BioStar, Zirkon BioStar Colour

- 12.4 Mobility in soil:** No information available.
- 12.5 Results of PBT and vPvB assessment:** The substances in the mixture do not meet the PBT/vPvB criteria according to REACH, annex XIII.
- 12.6 Other adverse effects:** No information available.

SECTION 13: Disposal considerations

- 13.1 Waste treatment methods**
- Advice on disposal:** Dispose of waste according to applicable legislation.
- Contaminated packaging:** Dispose of waste according to applicable legislation. Non-contaminated packages may be recycled.

SECTION 14: Transport information**Land transport (ADR/RID)**

- 14.1 UN number:** No dangerous good in sense of this transport regulation.
- 14.2 UN proper shipping name:** No dangerous good in sense of this transport regulation.
- 14.3 Transport hazard class(es):** No dangerous good in sense of this transport regulation.
- 14.4 Packing group:** No dangerous good in sense of this transport regulation.

Inland waterways transport (ADN)

- 14.1 UN number:** No dangerous good in sense of this transport regulation.
- 14.2 UN proper shipping name:** No dangerous good in sense of this transport regulation.
- 14.3 Transport hazard class(es):** No dangerous good in sense of this transport regulation.
- 14.4 Packing group:** No dangerous good in sense of this transport regulation.

Marine transport (IMDG)

- 14.1 UN number:** No dangerous good in sense of this transport regulation.
- 14.2 UN proper shipping name:** No dangerous good in sense of this transport regulation.
- 14.3 Transport hazard class(es):** No dangerous good in sense of this transport regulation.
- 14.4 Packing group:** No dangerous good in sense of this transport regulation.

Air transport (ICAO)

- 14.1 UN number:** No dangerous good in sense of this transport regulation.
- 14.2 UN proper shipping name:** No dangerous good in sense of this transport regulation.
- 14.3 Transport hazard class(es):** No dangerous good in sense of this transport regulation.
- 14.4 Packing group:** No dangerous good in sense of this transport regulation.

- 14.5 Environmental hazards**
ENVIRONMENTALLY HAZARDOUS: No
- 14.6 Special precautions for user:** No dangerous good in sense of this transport regulation.
- 14.7 Transport in bulk according to Annex II of Marpol and the IBC Code:** No dangerous good in sense of this transport regulation.

SECTION 15: Regulatory information

- 15.1 Safety, health and environmental regulations/legislation specific for the substance or mixture**
- National regulatory information**
Water contaminating class (D): - - not water contaminating
- 15.2 Chemical safety assessment:** For the following substances of this mixture a chemical safety assessment has been carried out:
Zirconia.

Zirkon BioStar, Zirkon BioStar Colour**SECTION 16: Other information**

Abbreviations and acronyms

ADR:	Accord européen sur le transport des marchandises dangereuses par Route (European Agreement concerning the International Carriage of Dangerous Goods by Road)
RID:	Règlement international conernat le transport des marchandises dangereuses par chemin de fer (Regulations Concerning the International Transport of Dangerous Goods by Rail)
IMDG:	International Maritime Code for Dangerous Goods
IATA:	International Air Transport Association
IATA-DGR:	Dangerous Goods Regulations by the "International Air Transport Association" (IATA)
ICAO:	International Civil Aviation Organization
ICAO-TI:	Technical Instructions by the "International Civil Aviation Organization" (ICAO)
CAS:	Chemical Abstracts Service (division of the American Chemical Society)
GHS:	Globally Harmonized System of Classification and Labelling of Chemicals
CLP:	Regulation on Classification, Labelling and Packaging of Substances and Mixtures,
LC50:	Lethal concentration, 50 percent
LD50:	Lethal dose, 50 percent
EC50:	Effective concentration, 50 percent
DNEL:	Derived No Effect Level
PNEC:	Predicted No Effect Concentration
PBT:	Persistent, Bioaccumulative and Toxic
vPvB:	very Persistent and very Bioaccumulative

Further Information

The above information describes exclusively the safety requirements of the product and is based on our present-day knowledge. The information is intended to give you advice about the safe handling of the product named in this safety data sheet, for storage, processing, transport and disposal. The information cannot be transferred to other products. In the case of mixing the product with other products or in the case of processing, the information on this safety data sheet is not necessarily valid for the new made-up material.

(The data for the hazardous ingredients were taken respectively from the last version of the sub-contractor's safety data sheet.)